IDEIA™ Norovirus

An enzyme immunoassay for the detection of Norovirus in human faecal specimens.

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IDEIA™

1 INTENDED USE

The IDEIA™ Norovirus test is a qualitative enzyme immunoassay for the detection of genogroup 1 and genogroup 2 Norovirus in human faeces.

2 SUMMARY

Norovirus, initially known as small round structured viruses^{1,2} (SRSV), are single stranded RNA viruses classified in the family Caliciviridae³.

Initially four antigenic types of SRSV were recognised⁴ but more recently, using molecular techniques, three genogroups have been identified within the genus Norovirus⁵. Genogroup 1 and Genogroup 2 are associated with human infections whilst Genogroup 3 is associated with bovine and porcine infection⁶.

Noroviruses are a major cause of viral gastroenteritis with high attack rates in children and adults^{7,8}. Large outbreaks can occur in hospitals⁹, cruise ships¹⁰, schools¹¹, and residential homes¹². Several outbreaks have arisen from the ingestion of contaminated water and food, especially shell fish. Outbreaks particularly of Norovirus associated with Winter Vomiting Disease can occur in the winter months in the western hemisphere but sporadic cases, and community and food borne outbreaks can occur throughout the year¹³.

Outbreaks of Norovirus are reported to be more common than bacterial gastroenteritis¹⁴ and can have a significant public health impact. Infections in infants, elderly patients and immunocompromised patients can give rise to severe dehydration. Outbreaks in hospitals can prolong hospitalisation and treatment of infected children.

The laboratory diagnosis of Norovirus infection plays an important role in patient and public health management and provides a differential diagnosis from other causes of gastroenteritis. At present Norovirus do not grow readily in cell culture systems and cannot be isolated from clinical specimens. Electron Microscopy has traditionally been used as an aid to diagnosis but is mainly confined to reference centres because of the high expertise and expense required. More recently a range of in-house molecular amplification methods have been described¹⁵ but these methods are not well standardised and performance is influenced by inhibitors present in specimens and by the primer sets used. In addition these tests are not convenient or cost effective for routine testing.

Recently genogroup and genotype specific antibodies^{16,17} have been generated and shown to be applicable for use in rapid enzyme immunoassays¹⁸. The IDEIA™ Norovirus test utilises genogroup 1 and genogroup 2 specific monoclonal and polyclonal antibodies in a solid phase immunoassay for the rapid detection of Norovirus genotypes.




3 PRINCIPLE OF THE TEST

The IDEIA™ Norovirus test utilises a combination of both genogroup 1 and genogroup 2 specific monoclonal and polyclonal antibodies in a solid-phase sandwich enzyme immunoassay to detect genogroup 1 and genogroup 2 Norovirus. Microwells are coated with a mixture of genogroup 1 and genogroup 2 Norovirus specific monoclonal antibodies. Faecal suspension or control specimen is added to the microwell and incubated simultaneously with a mixture of Norovirus specific genogroup 1 and genogroup 2 polyclonal and monoclonal antibodies conjugated to horseradish peroxidase.

Norovirus antigen present in the specimen is captured between antibodies on the solid phase and the enzyme conjugated antibodies. After 60 minutes incubation at room temperature, the microwells are washed with working strength wash buffer to remove excess specimen and any unbound enzyme labelled antibodies. A chromogen is added to the microwells and incubated for 30 minutes at room temperature. The presence of specifically bound enzyme labelled antibodies in the microwells results in a colour change which is stopped by the addition of acid. Colour intensity significantly above background levels is indicative of the presence of Norovirus antigen in the specimen or control. The IDEIA™ Norovirus test will detect both genogroup 1 and genogroup 2 strains, but will not differentiate between them.

4 DEFINITIONS

The following symbols have been used throughout the product information.

	Product code and catalogue number
	Consult the instructions for use
	Contains sufficient for 'N' tests
	Manufactured by
	<i>In vitro</i> diagnostic medical device
	Use by
	Batch Code
	Storage temperature limitations
	Add water

5 REAGENTS PROVIDED

 96 - Each kit contains sufficient materials for 96 direct specimens.  - The shelf life of the kit is as indicated on the outer box label.

5.1 IDEIA™ NOROVIRUS CONTENTS

	One Instructions for Use Booklet.
MICROTITRATION PLATE	1 x 96 well Microtitration Plate of 12, 8 microwell break-apart strips coated with a mixture of Norovirus genogroup 1 and 2 specific monoclonal antibodies. The plate is presented in a resealable foil pouch containing a desiccant card for subsequent storage of unused microwells.
One bottle of each of the following unless indicated otherwise:	
SAMPLE DILUENT	110mL Sample Diluent. Tris buffered saline containing antimicrobial agent and red dye.
CONTROL +	4 x Lyophilized Positive Control: of baculovirus expressed viral-like particles of Norovirus.

CONJUGATE	13mL Conjugate: peroxidase labelled polyclonal antibodies (rabbit) and monoclonal antibodies (mouse) – against Norovirus genogroup 1 and 2 in a buffered protein solution containing antimicrobial agent and blue dye.
CONTROL DILUENT	10mL Control Diluent: phosphate buffered saline containing antimicrobial agents and green dye.
WASH BUFFER (x25)	2 x 50mL Wash Buffer Concentrate (x25): tris buffered solution containing detergent and antimicrobial agent.
SUBSTRATE TMB	13mL Substrate: stabilized peroxide 3,3',5,5'-tetramethylbenzidine (TMB) in a dilute buffer solution. TMB has been reported to be non-carcinogenic. However, personal protective equipment is recommended to avoid direct exposure.
STOP SOLUTION	12mL Stop Solution: 0.46mol/L sulphuric acid.

5.2 PREPARATION, STORAGE AND RE-USE OF KIT COMPONENTS

The IDEIA™ Norovirus kit format allows for testing of up to 96 specimens. In order to ensure optimal kit performance it is important that unused kit components are stored according to the following instructions:

5.2.1 Antibody Coated Microwells - **MICROTITRATION PLATE**

Open the plate pouch by cutting along the seal. Break off the required number of microwells and relocate them into the frame. Return all unused microwells and strips to the resealable pouch with the desiccant card. Carefully reseal the pouch and store at 2-8°C. Microwells may be used for up to 12 weeks after initial opening, provided they are stored in this manner.

5.2.2 Sample Diluent - **SAMPLE DILUENT**

Ready to use. Store unused Sample Diluent at 2-8°C.

NOTE: Sample Diluent is required for use as a negative control and sufficient should be retained for use with each batch of specimens tested. When used as negative control sample diluent should be aliquoted into a vial identical to those used for preparation of specimen dilutions, to ensure that vials contain no interfering substances.

5.2.3 Positive Control - **CONTROL**

Reconstitute the required number of the lyophilised Positive Control with 1mL of Control Diluent. Mix the contents by inversion. Store unused positive control at 2-8°C and use within 1 month. For long term storage, store at -20°C in appropriate aliquots.

5.2.4 Conjugate - **CONJUGATE**

Ready to use. Store unused Conjugate at 2-8°C.

5.2.5 Control Diluent - **CONTROL DILUENT**

Ready to use. Store unused Control Diluent at 2-8°C.

5.2.6 Wash Buffer - [WASH BUFFER (X25)]

Provided x25 concentrate. Prepare working strength Wash Buffer by adding 1 part of Wash Buffer concentrate to 24 parts of fresh deionised or distilled water. **Prepare working strength Wash Buffer on the day of use.** Store unused concentrate at 2-8°C.

Do not store unused working strength Wash Buffer for subsequent use (see section 8.2.10).

5.2.7 Substrate - [SUBSTRATE TMB]

Ready to use. Store unused Substrate at 2-8°C, protect from light.

5.2.8 Stop Solution - [STOP SOLUTION]

Ready to use. Store unused Stop Solution at 2-8°C.

6 ADDITIONAL REAGENTS

6.1 REAGENTS

Fresh deionised or distilled water for preparation of working strength Wash Buffer.

6.2 ACCESSORIES

The following products are intended for use in conjunction with IDEIA™ Norovirus. Contact your local Oxoid subsidiary or distributor for further information.

IDEIA™ Gastroenteritis Sample Diluent 110mL (Code No. S6061).

IDEIA™ Gastroenteritis Wash Buffer Concentrate 50mL (Code No. S6126).

7 EQUIPMENT

The following equipment is required:

Containers suitable for the collection of faecal specimens.

Clean screw-capped disposable containers (minimum 3mL capacity) for preparation of faecal suspensions.

Clean vial for Negative Control.

Clean absorbent paper (onto which microwells can be tapped dry).

Precision micropipettes and disposable tips to deliver 100µL and 1000µL volumes (optional).

Waste discard container with suitable fresh disinfectant.

Automated plate washer (optional) or suitable equipment for washing 8 microwell strips (section 10.2.4).

NOTE: if washing less than 8 test microwells in a strip using an automated washer with an 8 microwell head, it is important to completely fill the strips with blank microwells.

Spectrophotometer or EIA plate reader able to read a 96 microwell plate of 8 microwell strips at an absorbance of 450nm with a reference at 620-650nm. (Optional, Section 10.3, Reading the Test Results).

Microplate mixer.

Application Notes for use on open automated systems are available for this assay. Contact your local Oxoid subsidiary or distributor.

8 PRECAUTIONS

 - For *in vitro* diagnostic use. Anyone performing an assay with this product must be trained in its use and must be experienced in laboratory procedures.

8.1 SAFETY PRECAUTIONS

8.1.1 Stop Solution contains sulphuric acid (0.46mol/L). Avoid contact with skin and mucous membranes. If stop solution comes into contact with these sites wash immediately with water.

8.1.2 Do not eat, drink, smoke, store or prepare foods, or apply cosmetics within the designated work area.

8.1.3 Do not pipette materials by mouth.

8.1.4 Wear disposable gloves when handling clinical specimens and reagents. Always wash hands after working with potentially infectious material.

8.1.5 Specimens, positive controls and all materials coming into contact with them should be handled and disposed of as though potentially infectious.

8.1.6 Dispose of all clinical specimens in accordance with local legislation.

8.1.7 IDEIA™ Norovirus reagents contain a proprietary anti-microbial agent which presents no hazard to the user if normal laboratory safety precautions are followed.

8.1.8 Safety data sheet available for professional user on request.

8.2 TECHNICAL PRECAUTIONS

8.2.1 Components must not be used after the expiry date printed on the labels. Do not mix or interchange different batches/lots of reagents.

8.2.2 The reagents are provided at fixed concentrations. Test performance will be affected if the reagents are modified or stored under conditions other than those detailed in Section 5.2.

8.2.3 Avoid contamination of reagents.

- 8.2.4** Use separate disposable pipettes or pipette tips for each specimen, control or reagent in order to avoid cross-contamination of either specimen or controls which could cause erroneous results.
- 8.2.5** When adding specimens or enzyme labelled antibody to the microwells, avoid allowing the pipette tip from coming into contact with microwell walls.
- 8.2.6** Avoid contamination with metal ions and oxidising agents.
- 8.2.7** Do not use substrate showing a blue colour prior to its addition to the microwells.
- 8.2.8** Protect substrate from light.
- 8.2.9** Microwells cannot be re-used.
- 8.2.10** Do not store unused working strength Wash Buffer for subsequent use. When not in use wash buffer reservoirs should be rinsed in deionised or distilled water and left to dry.
- 8.2.11** Manual or automated washing equipment must be free of microbial contamination, be correctly calibrated and maintained according to the manufacturer's instructions.

9 COLLECTION AND PREPARATION OF SPECIMENS

9.1 SPECIMEN COLLECTION

Faecal specimens should be collected as soon as possible following the onset of symptoms.

Peak excretion of Norovirus in faeces from patients with gastroenteritis is reported to occur between 25-72 hours after the onset of symptoms.

Faecal specimens for direct testing should be collected into containers that do not contain media, preservatives, animal sera, metal ions, oxidising agents or detergents, as all of these additives may interfere with the IDEIA™ Norovirus test.

Specimens may be stored for 3 days at 2-8°C prior to testing. For long-term storage of faecal specimens, store at -20°C.

9.2 PREPARATION OF FAECAL SPECIMENS

Add 1mL of Sample Diluent to a suitable labelled container and use to prepare a 10% suspension or dilution of faecal specimen by addition of approximately 0.1g of solid faeces (small pea-sized portion) or approximately 100µL of liquid faeces. Mix thoroughly and leave to settle for 10 minutes prior to testing.

Alternatively, after preparation of samples centrifuge the suspension for 5 minutes at ≥5000 rpm.

Specimens suspended/diluted in IDEIA™ Norovirus Sample Diluent may be stored at 2-8°C for up to 3 days prior to testing.

NOTE: Faecal specimens prepared in Amplified IDEIA™ Astrovirus, IDEIA™ Rotavirus and IDEIA™ Adenovirus Sample Diluent can also be tested in the IDEIA™ Norovirus test. Alternative sample diluents have not been validated for use.

10 TEST PROCEDURE

PLEASE REFER TO SECTION 8.2, TECHNICAL PRECAUTIONS, BEFORE PERFORMING TEST PROCEDURE.

Reagents and specimens should be brought to room temperature (15-30°C) before use.

10.1 PREPARATION OF CONTROLS

Negative Control

Add 1mL of Sample Diluent to a vial identical to those used for specimen dilution.

Positive Control

The Positive Control must be reconstituted with 1mL of Control Diluent and mixed before use (see Section 5.2.3).

At least one IDEIA™ Norovirus Positive and Negative Control must be included with each batch of specimens tested (see Section 10.2.1).

10.2 ASSAY PROCEDURE

When using reagent dropper bottles, hold the bottles vertically with the nozzle approximately 5mm above the microwell. Squeeze the bottle gently and ensure that the drops fall freely to the microwells without touching the sides. Avoid contamination of all the dropper bottle nozzles. Specimens, controls and conjugate should be added using a micropipette.

10.2.1 Specimen Addition

Locate the required number of microwells into the microwell holder.

Add 100µL of each diluted faecal specimen/stool supernatant, Negative Control or Positive Control to separate microwells. (At least one Negative Control and one Positive Control should be included in each batch of tests.)

Avoid disturbing or sampling any sediment in the sample vial. When dispensing specimens, hold the transfer pipette vertically with the tip immediately above the microwell and ensure that the specimen is delivered without touching the sides of the microwell. Take care not to cross-contaminate specimen and control microwells as this may cause erroneous results.

10.2.2 Conjugate Addition

NOTE: If multiple strips of 8 microwells are used it is recommended to use an 8-channel pipette.

After addition of all specimens and controls to the microwells.

Add 100µL of Conjugate to all microwells using a micropipette and mix gently.

10.2.3 First Incubation

Incubate all microwells 20-30°C for **60 ± 5 minutes**.

10.2.4 Washing the Microwells

The microwells should be washed using freshly prepared working strength Wash Buffer (see Section 5.2.6).

The washing technique is critical to the test performance (see Section 8.2.11) and should be carried out so as to ensure complete filling (with a minimum of **350µL** of working strength Wash Buffer) and emptying of the microwells.

Five wash cycles are essential, by either automated or manual washing techniques.

Manual Washing

If washing microwells manually, aspirate or shake out the contents of the microwells and using freshly prepared Wash Buffer, ensure complete filling and emptying of microwells. Between each wash step remove all remaining Wash Buffer by tapping the inverted microwells on to clean absorbent paper. Manual washing efficiency is further ensured if the Wash Buffer is delivered at an angle so as to produce a vortex in the microwells. After the final wash, the plate should be inverted and tapped on absorbent paper to remove the last traces of Wash Buffer.

Automated Washing

Automated washers should be programmed to complete 5 wash cycles. Washers must be correctly calibrated to ensure complete filling and emptying of microwells between each wash. After the final wash, the plate should be inverted and tapped on absorbent paper to remove the last traces of Wash Buffer.

10.2.5 Substrate Addition

Add 2 drops (or 100µL) of Substrate to each microwell.

10.2.6 Second Incubation

Incubate at 20-30°C for **30 minutes**.

Microwells can be **read visually immediately** after the second incubation (see Section 10.3.1).

10.2.7 Stopping the reaction

Add 2 drops (or 100µL) of Stop Solution to each microwell.

Ensure thorough mixing of the microwells before reading the results. The coloured product is stable for up to 30 minutes after addition of stop solution.

10.3 READING THE TEST RESULTS

10.3.1 Visual reading


The microwells may be assessed visually immediately after the second incubation. It is recommended that any microwells in which the colour intensity is difficult to interpret when compared to the Negative Control are also read photometrically after addition of Stop Solution (see Section 10.3.2).

10.3.2 Photometric reading

The microwells should be read photometrically within 30 minutes after addition of the Stop Solution. Mix the contents of the microwells and read the absorbance of each microwell using a suitable spectrophotometer or EIA plate reader set at 450nm. Ensure that the bottoms of the microwells are clean before reading and check that no foreign matter is present in the microwells. The reader should be blanked on air (ie with no plate in the carriage) before the plate is scanned.

Alternatively, if the spectrophotometer allows for the use of a reference wavelength (at 620 to 650nm), dual wavelength reading should be performed as this eliminates any potential interference caused by aberrations, such as dirt or marks, on the optical surface of the microwells.

10.4 SUMMARY OF THE IDEIA™ NOROVIRUS ASSAY PROCEDURE

11 INTERPRETATION OF TEST RESULTS

11.1 NEGATIVE CONTROL

As detailed in Section 10.1 (Preparation of Controls), at least one Negative Control must be included in each assay run.

Photometric determination

The Negative Control value, or mean of the Negative Control values, should be less than 0.150 absorbance units.

11.2 CUT-OFF VALUE

Calculation of cut-off value

The cut-off value is calculated by adding 0.10 absorbance units to the negative control value or mean value when more than one Negative Control is included.

Visual determination

The Negative Control microwell should have no visible colour. If this is not the case the test should be read photometrically or repeated.

11.3 POSITIVE CONTROL

As detailed in Section 10.1 (Preparation of Controls), at least one Positive Control must be included in each assay run.

Visual determination

The Positive Control microwell should show a distinct blue colour clearly distinguishable from the Negative Control. If this is not the case test results should not be determined visually. Results should be read photometrically or the test repeated.

Photometric determination

The positive control must have a value of greater than 0.5Au.

11.4 SPECIMENS

11.4.1 Visual determination

Any specimen showing a blue colour more than the Negative Control value should be interpreted as positive, provided that the Negative Control is colourless. Any specimens with no visible colour should be considered negative. If the Negative Control shows any visible blue colour results must be read photometrically after addition of Stopping Solution or the test repeated.

If, after the addition of Substrate, the content of a microwell turns dark blue and forms a blue-black precipitate, the specimen should be interpreted as positive.

11.4.2 Photometric determination

Any specimen with an absorbance value greater than the cut-off value is positive (see Section 11.1). Any specimen with an absorbance value less than the cut-off should be considered negative. A result within 0.010 absorbance unit of the cut-off value should be considered equivocal, and the test repeated or the patient resampled.

12 PERFORMANCE LIMITATIONS

12.1 A negative result does not exclude the possibility of norovirus infection in the patient. Failure to detect norovirus may be a result of factors such as collection of specimen at an improper time in the disease when too few virions are present, improper sampling and/or handling of the specimen.

12.2 IDEIA™ Norovirus test detects genogroup specific viral proteins present in human genotypes of genogroup 1 and genogroup 2 Norovirus. The test cannot be used to differentiate between serotypes of Norovirus.

12.3 The reagents are provided at fixed working concentrations. Test performance will be adversely affected if reagents are modified or stored under conditions other than those detailed in Section 5.2.

12.4 All positive results must be interpreted in conjunction with patient related clinical information. Test results should be interpreted in conjunction with information available from epidemiological studies, clinical assessment of the patient and other diagnostic procedures.

12.5 Meconium samples have not been validated for use with the IDEIA™ Norovirus.

12.6 The IDEIA™ Norovirus test has not been validated with all known strains of Norovirus and therefore failure to detect Norovirus may be due to antigenic diversity in circulating strains.

12.7 A positive result does not preclude the presence of other enteric pathogens. Whilst the relationship between Norovirus and gastroenteritis is well established, concurrent infection with other microbial pathogens is possible. Additional microbiological tests should be performed in parallel with the IDEIA™ Norovirus test in order to exclude other possible causes of illness.

13 EXPECTED VALUES

Positivity rates may vary according to the prevalence of Norovirus in different populations, genotype circulating, geographical location, method of specimen collection, handling, storage, transportation of specimens and the general health environment of the patient population under study.^{7,13}

Noroviruses are a common cause of community outbreaks of gastroenteritis particularly in hospitals or residential homes associated with >40% of outbreaks in the UK^{9,12}.

Noroviruses are frequently associated with cases of winter vomiting disease.

14 SPECIFIC PERFORMANCE CHARACTERISTICS

14.1 CLINICAL STUDIES

The Updated IDEIA™ Norovirus test was evaluated against the previous version of the test (K6043). Studies were conducted on faecal specimens collected from outbreaks of gastroenteritis suspected of being caused by Norovirus from across the UK. All specimens had previously been tested by PCR. Results of the two IDEIA™ Norovirus assays were compared with the known PCR status of the specimens.

14.2 CLINICAL PERFORMANCE

A total of 120 specimens were tested comprising 83 specimens with a positive PCR result for Norovirus and 37 specimens with a negative PCR result for Norovirus.

The updated IDEIA™ Norovirus assay had improved sensitivity and negative predictive value (NPV) relative to PCR, and equivalent specificity and positive predictive value (PPV) relative to PCR when compared with the previous version of the kit.

Table 14.2 Performance relative to PCR of updated IDEIA™ Norovirus (K6044) and the previous version IDEIA™ Norovirus (K6043).

	Updated IDEIA™ Norovirus K6044	Previous version IDEIA™ Norovirus K6043
Sensitivity relative to PCR	72.8% (59/81)	55.4% (46/83)
Specificity relative to PCR	100% (37/37)	100% (37/37)
PPV	100% (59/59)	100% (46/46)
NPV	62.7% (37/59)	51.4% (83/120)
Outbreaks detected (No. detected by PCR)	16 (21)	12 (21)

PPV = Number of true positive samples correctly identified by the test

NPV = Number of true negative samples correctly identified by the test

14.3 CROSS REACTIVITY

The following micro-organisms were tested and shown to be negative in the IDEIA™ Norovirus test. Cross-reactivity tests were performed either on clinical specimens for which the microbial status had been determined, or on laboratory cultures of known organisms, containing approximately 10^7 - 10^8 viable organisms/mL. The source of micro-organisms is referenced in the key below:

Viruses

Adenovirus^a

Astrovirus^a

Rotavirus^a

Peptococcus sp

Peptostreptococcus sp

Proteus sp

Pseudomonas sp

Salmonella enteritidis^a

Salmonella typhimurium

Serratia sp

Shigella flexner^a

Shigella sonnei

Staphylococcus aureus

Staphylococcus epidermidis

Streptococcus group A

Veillonella sp

Vibrio cholerae^a

Vibrio parahaemolyticus^a

Bacteria

Acinetobacter sp

Aeromonas sp

Bacillus sp

Campylobacter coli^a

Campylobacter jejuni^a

Citrobacter sp

Clostridium perfringens^a

Clostridium difficile^a

Enterobacter sp

Escherichia coli

Enterococcus faecalis

Haemophilus influenzae

Helicobacter pylori

Klebsiella sp

Lactobacillus sp

Listeria sp

Other micro-organisms

Candida albicans^a

Key:

^aMicro-organisms present and tested in faeces

^bMicro-organisms tested both in faeces and in broth culture

All other micro-organisms were tested in broth culture

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1. **Jiang, X., D.Y. Graham, K. Wang, and M.K. Estes. (1990)**
Norwalk virus genome cloning and characterisation
Science **250**:1580-1583.
2. **Jiang, X.M. Wang, K. Wang, and M.K. Estes (1993)**
Sequence and genomic organisation of Norwalk virus.
Virology **195**:51-61.
3. **Green, K. Y., T. Ando, M.S. Balayan, I.N. Clarke, M.K. Estes, D.O. Matson, S. N Nakata, J.D. Neill, M.J. Studdert, and H.J. Theil**
Caliciviridae, In M.H. van Regenmortel, C.M. Fauquet, D.H.L. Bishop, E.B. Carsten, M.K. Estes, S.M. Lemon, J. Maniloff, M.H. Mayo, D.J. McGeoch, C.R. Pringle, and R.B. Wickner (ed.), *Virus taxonomy*. Seventh report of the International Committee on Taxonomy of Viruses, in press. Academic Press, Orlando, Fla.
4. **Lewis, D. (1991)**
Norwalk agent and other small round structured viruses in the U.K.
J. Infect. **23**:220-222.
5. **Green S.M., K.E. Dingle, P.R. Lambden, E.O. Caul, C.R. Ashley, and I.N. Clarke (1994)**
Human enteric Caliciviridae; a new prevalent SRSV group defined by RNA-dependent RNA polymerase and capsid diversity.
J. Gen. Virol. **75**:1883-1888.
6. **Ando T, Noel JS, Fankhauser (2000)**
RL. Genetic classification of "Norwalk-like viruses."
J Infect Dis **181**(suppl 2):S336—48.
7. **Blacklow, N.R. and H.B. Greenberg. (1990)**
Viral gastroenteritis.
N. Engl. J. Med. **325**:252-264.
8. **E. Owen Caul (1996)**
Viral gastroenteritis: Small round structured viruses, caliciviruses and astroviruses. Part 1. The clinical and diagnostic perspective.
J Clin Pathol **49**:874-880.
9. **Chadwick, P.R., and R. McCann, (1994)**
Transmission of a small round structured virus by vomiting during a hospital outbreak of gastroenteritis.
J. Hosp. Infect. **26**:251-259.
10. **McEvoy, M., W. Balke, D. Brown, J. Green, and R. Cartwright. (1996)**
An outbreak of viral gastroenteritis on a cruise ship.
Commun. Dis. Rep. CDR Rev. **6**:R188-R192.
11. **Kobayashi, S., T. Morishita, T. Yamashita, K. Sakae, O. Nishio, T. Miyake, Y. Ishihara, and S. Isomura. (1991)**
A large outbreak of gastroenteritis associated with a small round structured virus among schoolchildren and teachers in Japan.
Epidemiol. Infect. **107**:81-86.

12. **Jiang, X., E. Turf, J.Hu, E. Barrett, X.M. Dal, S. Monroe, C. Humphrey, L.K. Pickering, and D.O. Matson. (1996)**
Outbreaks of gastroenteritis in elderly nursing homes and retirement facilities associated with human caliciviruses.
J. Med. Virol. **50**:335-341.
13. **E. Owen Caul (1996)**
Viral gastroenteritis: Small round structured viruses, calciviruses and astroviruses. Part II. The epidemiological perspective.
J Clin Pathol.**49**:959-964.
14. **H.S Evans, P. Madden, C. Douglas, G.K. Adak, S.J. O'Brien, T. Djuretic, P.G. Wall, R. Stanwell-Smith. (1998)**
General outbreaks of infectious intestinal disease in England and Wales: 1995 and 1996.
Commun Dis Public Health **1**:165-71.
15. **Green J. Hale A.D, Brown DWG. (1995)**
Recent developments in the detection and characterisation of small round structured viruses.
PHLS Microbiology Digest:**12**:219-22.
16. **James P. Brinkler, Neil R. Blacklow, Mary K. Estes, Christine L. Moe, Kellogg J. Schwab, and John E. Herrmann. (1998)**
Detection of Norwalk Virus and other genogroup 1 human Caliciviruses by a Monoclonal Antibody, Recombinant-Antigen-Based Immunoglobulin M. Capture Enzyme Immunoassay.
Journal of Clinical, Microbiology.p.1064-1069.
17. **Antony D. Hale, Tomoyuki N. Tanaka, Noritoshi Kitamoto, Max Ciarlet, Xi Jiang, Naokazu Takeda, David W.G. Brown, and M. K. Estes. (2000)**
Identification of an Epitope Common to Genogroup 1 "Norwalk-Like Viruses".
Journal of Clinical Microbiology. p.1656-1660.
18. **Vipond IB, Pelosi E, Williams J, Ashley CR, Lambden PR, Clarke IN, Caul EO. (2000)**
A diagnostic EIA for detection of the prevalent SRSV strain in United Kingdom outbreaks of gastroenteritis.

TECHNICAL ADVICE AND CUSTOMER SERVICE
CONSEIL TECHNIQUE ET SERVICE CLIENTELE
TECHNISCHE BERATUNG UND KUNDENDIENST

For all enquiries please contact your local Oxoid.
subsidiary or distributor.

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