
**REF****M80 Rotascreen®****IVD****INTENDED USE**

Rotascreen® is a rapid and simple latex agglutination test for the detection of Rotavirus antigen in faecal specimens. The test is intended for professional use only.

**PRINCIPLE OF THE TEST**

Latex particles in the Rotascreen® Test Reagent are coated with rabbit antibodies raised against a pool of different Rotavirus isolates, including human.

When a faecal extract is mixed with the Test Reagent any Rotavirus antigens present will react with the sensitising antibodies, resulting in visible agglutination of the latex particles.

A Control Reagent, latex particles coated with normal rabbit globulin, is included to identify non-specific reactions which may occur with some faecal specimens.

**CONT**

<b>REAG TEST</b>	M80a	2.5mL
------------------	------	-------

Test reagent - Contains (rabbit) Rotavirus antibody-sensitised latex particles in buffer, with stabiliser, and 0.099% sodium azide as preservative.

<b>REAG CONTROL</b>	M80b	2.5mL
---------------------	------	-------

Control Reagent - Contains (rabbit) normal globulin-sensitised latex particles in buffer, with stabiliser, and 0.099% sodium azide as preservative.

<b>CONTROL +</b>	M80c	1.0mL
------------------	------	-------

Positive Control - Contains inactivated bovine Rotavirus antigens in buffered cell culture medium, with antibiotics, and 0.099% sodium azide as preservative.

<b>SAMP DIL</b>	M80d	50mL
-----------------	------	------

Sample Diluent - Contains working strength buffer, pH 7.2, and 0.099% sodium azide as preservative

Disposable Reaction Cards

Disposable Mixing Sticks

**ADDITIONAL MATERIALS REQUIRED BUT NOT PROVIDED**

Laboratory centrifuge

Stoppered or screw-capped tubes

**WARNINGS AND PRECAUTIONS**

1. Rotascreen® is for *in vitro* diagnostic use only.
2. Do not use reagents after the kit expiration date.
3. Do not interchange reagents from different kit lots.
4. The test should only be performed in accordance with the instructions supplied with the kit.
5. Do not pipette specimens or reagents by mouth.
6. All clinical specimens should be considered infectious and handled and disposed of according to accepted practices. Decontamination of infectious material can be achieved with sodium hypochlorite at a final concentration of 3% for 30minutes. Liquid waste containing acid must be neutralised before treatment.
7. The Positive Control is prepared from inactivated Rotavirus. However, it should be considered as potentially infectious and handled in the same way as a clinical specimen.

8. The reagents in this kit contain sodium azide as preservative which may react with lead or copper plumbing to form highly explosive azides. Upon disposal, flush with a large volume of water to prevent azide build-up.

**STORAGE**

Store all reagents at 2-8°C. **Do not freeze.** Under these conditions the reagents will be usable until the date printed on the outer carton label.

**INDICATIONS OF DETERIORATION**

Deterioration should be suspected if:-

- Clumping of the Test Reagent (M80a) or Control Reagent (M80b) is evident and cannot be removed by gentle mixing.
  - The Positive Control (M80c) or Sample Diluent (M80d) becomes cloudy or forms a sediment.
  - The Positive Control (M80c) fails to cause agglutination of the Test Reagent within the recommended reaction time.
- Reagents showing signs of deterioration should not be used.

**SPECIMEN STORAGE AND PREPARATION**

Faecal specimens should be tested as soon after collection as possible but, if necessary, they may be stored overnight at 2-8°C or for longer periods at -20°C or below. Contaminated samples may give erroneous results.

Prepare an approximate 10% suspension of the faecal specimen by transferring 0.1mL (0.1g) of sample into 1.0mL of **Sample Diluent** in a stoppered or screw-capped tube. Mix the contents well.

Allow to stand at Room Temperature for 1-2 minutes before proceeding with the test.

**USE OF THE POSITIVE CONTROL**

The Positive control should be tested regularly to ensure that the reagents are functioning correctly. The control is ready for use and should be tested in place of the specimen extract in the Test Procedure.

The Positive control should give agglutination of the Test Reagent with no agglutination of the Control

Reagent. Failure to give this agglutination pattern is evidence of reagent deterioration.

**TEST PROCEDURE - STANDARD METHOD USING CENTRIFUGE**

1. Process the specimen as described in Section 8.
2. Allow the reagents to reach room temperature.
3. Centrifuge specimen extract at approximately 1000g for 10 mins.
4. Pipette 50µL of clear supernatant onto each of two wells on the reaction card.
5. Add 1 drop of well-mixed Test Reagent (M80a) to one well and 1 drop of well-mixed Control Reagent (M80b) to the other.
6. Mix the contents of each well using a separate mixing stick for each sample, covering the entire area of the well.
7. Gently rock the card and observe the agglutination for up to two minutes.

**TEST PROCEDURE - FILTER METHOD USING MICROGEN BIOPRODUCTS FILTER PACK M802**

1. Process the specimen as described in Section 8.
2. Allow the reagents to reach room temperature.
3. Remove the stopper and fit integral filter/dropper unit.

4. Holding the whole assembly vertically dispense 1 drop of clear filtrate onto each of two wells on the Test slide.
5. Add 1 drop of well-mixed Test Reagent (M80a) to one well and 1 drop of well-mixed Control Reagent (M80b) to the other.
6. Mix the contents of each well using a separate mixing stick for each sample, covering the entire area of the well.
7. Gently rock the card and observe the agglutination for up to two minutes.

#### INTERPRETATION OF RESULTS

- a) A **positive** result is indicated by agglutination of the Test Reagent (M80a) with no agglutination of the Control Reagent (M80b).
- b) The result is **negative** if no agglutination of either the Test Reagent (M80a) or the Control Reagent (M80b) is observed within the 2 minute test period.

Note: agglutination of the Control Reagent (M80b) is evidence of a non-specific reaction and means that the specimen is unsuitable for testing by this method.

#### LIMITATIONS OF THE PROCEDURE

Rotascreen® results should be evaluated in the light of all other available clinical and laboratory information. A positive Rotascreen® test does not preclude the possibility of other microbial co-infections. Rotascreen® is intended as an acute-phase test. Faecal samples collected after the acute phase may contain antigen concentrations below the threshold of reagent sensitivity.

#### PERFORMANCE CHARACTERISTICS

Rotascreen® has been independently evaluated in clinical use at three leading microbiology laboratories in the United Kingdom. Their results can be summarised as follows:

##### Laboratory 1

69 stool specimens were tested by Rotascreen®, in parallel with a respected commercial enzyme immunoassay (EIA) system:

**Table 1 Comparison of Rotascreen® and a commercial EIA system**

	Rotascreen®	Rotascreen®	Total
	+	-	
EIA +	36	1*	37
EIA -	0	32	32
Total	36	33	69

\*EM negative

Sensitivity	97%
Specificity	100%
Predictive value positive	100%
Predictive value negative	97%

##### Laboratory 2

398 stool specimens were tested for Rotavirus by Rotascreen®, electron microscopy, polyacrylamide gel electrophoresis (PAGE) RNA banding, and the WHO reference EIA system.

**Table 2 Comparison of Rotascreen®, EIA, Electron Microscopy and PAGE electrophoresis**

Rotascreen®	EIA	EM	PAGE	No. of Samples	% of total samples
+	+	+	+	161	40.25
+	+	+	-	18	4.5
-	+	+	+	8	2.0
+	+	-	+	22	5.5
-	+	-	+	3	0.75
-	+	+	-	7	1.75
+	+	-	-	62	15.5
-	+	-	-	22	5.5
-	-	-	-	95	23.75

If a true Rotavirus positive is defined as a specimen which is positive by two or more tests, performance of the different assays can be summarised as follows:

**Table 3 Summary of performance of Rotascreen®, in comparison with EIA, EM and PAGE**

	Rotascreen®	EIA	EM	PAGE
Sensitivity	93%	100%	69%	69%
Specificity	100%	81%	100%	100%
Predictive value positive	100%	93%	100%	100%
Predictive value negative	87%	100%	57%	57%

##### Laboratory 3

278 stool specimens were characterised by electron microscopy, EIA and PAGE RNA banding. True positives were considered to be those specimens positive by one or more of these tests. True negatives were those which were negative by all tests. The characterised sera were then tested by Rotascreen® in parallel with two other commercially available slide latex tests. Results can be summarised as follows:

**Table 4 Comparison of Rotascreen® and two commercially available slide latex tests**

	Rotascreen®	Latex 1	Latex 2
Sensitivity	91%	87%	76%
Specificity	97%	94%	86%
% of specimens which agglutinated both test and control latex reagents	1%	2%	14%
Predictive value positive*	98%	97%	99%
Predictive value negative*	95%	90%	92%

- \* Note: specimens giving agglutination of both Test and Control latex reagents were excluded from the calculation of Predictive Values.

##### Intra Batch Reproducibility

Intra batch reproducibility has been determined by running a panel of 5 positive and 5 negative samples on 10 separate occasions using four different operators. No differences in interpretation was noted on any of the test dates or across any of the operators.

**Interbatch sensitivity and specificity** is maintained by monitoring each batch against reference panels of samples and serial dilutions of a sensitivity standard (Rota positive sample).

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Microgen Bioproducts Limited