# Combi K-SeT

www.corisbio.com

IFU-5804/T1/04

## In vitro rapid diagnostic test for the detection of Rotavirus and Adenovirus in human stool sample

FOR IN VITRO USE  
FOR PROFESSIONAL USE ONLY

References: K-1504, 20 tests per kit, with collection set  
K-1204, 20 tests per kit, without collection set

### I. INTRODUCTION

Diarrhoea and gastro-enteritis in human beings can be caused by viruses (Rotavirus, Adenovirus, Astrovirus, calicivirus, etc), bacteria such as Salmonella and E. coli, and protozoa such as Cryptosporidium and Giardia. Viruses cause 45% of the diarrhoea in children under 1 year old and 40% of the diarrhoea in children under four years.

Rotavirus is the leading cause of gastro-enteritis in children under five years. Adenovirus's prevalence is 4-12%, which puts it in second place as a cause of viral enteritis in children under two years.

Rotavirus is transmitted by faecal-oral contact. After an incubation period of about three days it triggers fever, vomiting, and diarrhoea that can persist for up to ten days. It is responsible for 140 million cases of diarrhoea per year, with 870,000 deaths mainly due to dehydration in developing countries (WHO 1997). It is thus a major cause of mortality in Third World countries. However, it remains a severe infection even in the developed world. In the United States, 75 to 125 children are estimated to die from rotavirus infections each year. As it is highly contagious, it spreads very rapidly in paediatric populations, which are risk groups for such infections.

Adenovirus infection occurs by the faecal-oral route, but can also result from inhalation. The incubation period is from five to eight days and the symptoms of the stomach and intestinal inflammation are watery diarrhoea, vomiting, fever, and abdominal cramps.


The Adenoviruses are divided into six subgroups labelled A to F. Subgroup F is the most frequently involved in paediatric gastro-enteritis. The Combi K-SeT detects all Rotavirus and Adenovirus group viruses.

### II. PRINCIPLE OF THE TEST

This test is ready to use and is based on the homogeneous membrane system technology with colloidal gold particles. A nitrocellulose membrane is sensitized with antibodies directed against Rotavirus and Adenovirus (test lines). Test's specificity is due to two monoclonal antibodies directed against Rotavirus and specific proteins of human Adenovirus (Hexon antigen), respectively, that are conjugated to colloidal gold. These conjugates are insolubilized on a polyester membrane.

Faecal sample must be diluted in the dilution buffer that is supplied with the test. When 4 drops of the liquid phase of the faecal suspension come into contact with the strip, the solubilised conjugate migrates with the sample by passive diffusion and both conjugate and sample material come into contact with monoclonal antibody directed against specific Adenovirus proteins. If the sample contains Adenovirus, conjugate-Adenovirus complex will remain bound to the monoclonal antibody adsorbed onto the nitrocellulose and a red line will develop. Solution continues to migrate to encounter the anti-Rotavirus monoclonal that is adsorbed onto the nitrocellulose. If the sample contains Rotavirus, conjugate-Rotavirus complex will remain bound to the anti-Rotavirus monoclonal and a red line will develop.

Solution continues to migrate to encounter a third reagent (an anti-chicken IgY polyclonal antibody) that binds the migration control conjugate, thereby producing a red control line that confirms that the test is working properly. Result is visible within 10 minutes.

Manufacturer:

Coris BioConcept  
Science Park CREALYS  
B – 5032 GEMBLOUX  
BELGIUM  
Tel.: +32(0)81.719.917  
Fax: +32(0)81.719.919  
info@corisbio.com

Produced in BELGIUM

### III. REAGENTS AND MATERIALS

#### 1. Combi K-SeT (20)

- 20 sealed pouches containing one device.  
Each device contains one sensitized strip. Each device is individually packed in a pouch with a desiccant.

#### 2. Instruction for use (1)

#### 3. HC dilution buffer

Saline dilution buffered to pH 7.5 containing Tris, EDTA, NaN<sub>3</sub> (<0,1%), a detergent and blocking proteins.

- K-1204: 1 vial (15 mL)
- K-1504: 20 Faecal Sampling System (FSS) (2 mL) with a sampling screw

#### Materials to be ordered separately:

- Rotavirus positive control (Ref.: C-1081)
- Adenovirus positive control (Ref.: C-1082)
- Negative control (Ref.: CTR-1000)

### IV. SPECIAL PRECAUTIONS

- All operations linked to the use of the test must be performed in accordance with Good Laboratory Practices (GLP).
- All reagents are for *in vitro* diagnostic use only.
- Pouch must be opened with care
- Avoid touching nitrocellulose with your fingers.
- Wear gloves when handling samples.
- Never use reagents from another kit.
- Green lines indicate immunoreagents adsorption sites. Green color disappears during the test.
- Reagents' quality cannot be guaranteed beyond their shelf-life dates or if reagents are stored under inappropriate conditions.

### V. WASTE DISPOSAL

- Dispose of gloves, swabs, test tubes and used devices in accordance with GLP.
- Each user is responsible for the management of any waste produced, and must ensure that it is disposed of in accordance with applicable legislation

### VI. STORAGE

An unopened pouch may be kept at between 4 and 30°C and used until shelf-life date indicated on the packaging. Once pouch is opened, run test immediately.

Avoid freezing devices and buffer.

### VII. SPECIMEN HANDLING AND COLLECTION

The stool specimens must be tested as soon as possible after they are collected. If necessary, they may be stored at 2-8°C for 1 week or -20°C for longer periods of time.

Make sure that the specimens are not treated with solutions containing formaldehyde or its derivatives.

### VIII. PROCEDURE


#### Preparations of the test:

Allow kit components, in unopened packaging, and specimens to reach the room temperature (15-30°C) before performing a test.

Open pouch and remove the device. Once opened, run test immediately. Write patient's name or specimen number on device (one device per sample).

#### Specimen preparation procedure with Faecal Sampling System (K-1504):

1. Open the FSS tube and use the screw to collect stool (A). **The dilution ratio must be at most 4% w/v.** Take care to not take too much (B1) or too few specimens (B2). For liquid or semi-liquid samples, pour 80 µL of sample using a plastic pipette (not provided) into the FSS vial.
2. Insert the screw into the FSS and tighten the cap. Vortex the preparation to homogenize (C). The entire stool sample must be diluted into the solution.
3. Break off the point of the cap (D) and dispense 4 drops of diluted sample into the sample well of the device as illustrated below. To assure proper delivery, FSS vial must be held vertically (E).

Specimen preparation procedure (K-1204):

1. Add 0.5 mL or 15 drops of the dilution buffer solution into a tube.
2. Dip a loop containing the stool sample into the tube. **The dilution ratio must be at most 4% w/v.** For liquid samples, take 2 loops of 10 µL, for solid samples, take 1 loop.
3. Vortex the preparation to homogenize. The entire stool sample must be diluted into the solution.
4. Slowly dispense 100µL of diluted sample into the sample well of the device as illustrated below.

Let react for 10 minutes. Results must be observed in reading window. Positive results may be reported sooner as long as the test and control lines become visible.

**Result must be read on still wet strip**


**IX. INTERPRETING RESULTS**

Results are to be interpreted as follows:

**Negative test result:** a reddish-purple band appears across the central reading window at the Control line (C). No other band is present.

**Positive test result:** in addition to a reddish-purple band at the Control line (C), a visible reddish-purple band appears at the Rotavirus test line position (R) or Adenovirus test line (A). 3 lines (C-R-A) will appear in the case of an infection by both Rotavirus and Adenovirus. Intensity of the test line may vary according to the quantity of antigens found in the sample. Any reddish-purple test line (R or A), even weak, should be considered as a positive result.

**Invalid test result:** No detectable reddish-purple band at the Control line (C). The absence of a Control line indicates a failure in the test procedure, Repeat Invalid tests with a new test device.

Note: during the drying process, a very faint shadow may appear at the Test line position. It should not be regarded as a positive result.

**X. QUALITY CONTROL**

In accordance with Good Laboratory Practices, we recommend to check the test's performance regularly according to the laboratory's requirements. For the test, consider positive control as a liquid stool (see VIII.1).

**XI. PERFORMANCE (based on Combi-Strip kit)**

There is an excellent agreement (100%) between Combi K-Set kit and standard Combi-Strip kit.

**A. Detection Limit**

For Adenovirus, analytical sensitivity has been realised with a purified adenovirus-5 antigen and has been evaluated at about 3.9 x 10<sup>8</sup> vp/mL.

**B. Sensitivity - Specificity (Correlation):**

An evaluation has been carried out in which results obtained with Combi kit were compared with those of two marketed ELISA tests (for Rotavirus and for Adenovirus group).

Rotavirus and Adenovirus Combi kit's sensitivities and specificities were tested on 214 and 130 stool samples, respectively. Following results were obtained:

ELISA Rota	Positive	Negative	Total
<b>Coris BioConcept</b>			
Positive	105	0	105
Negative	1	108	109
<b>Total</b>	<b>106</b>	<b>108</b>	<b>214</b>

Sensitivity for Rotavirus: 99.1% (105/106)

Specificity for Rotavirus: 100% (108/108)

Reliability (Concordance): 99.5% (213/214)

ELISA Adenovirus	Positive	Negative	Total
<b>Coris BioConcept</b>			
Positive	12	0	12
Negative	0	118	118
<b>Total</b>	<b>12</b>	<b>118</b>	<b>130</b>

Sensitivity for Adenovirus: 100%

Specificity for Adenovirus: 100%

Reliability (Concordance): 100% (130/130)

**C. Accuracy**

To check intra-batch accuracy, same positive samples and a buffer solution have been processed 15 times on Combi kits of the same production batch in the same experimental conditions. All observed results were correct as expected. All 15 tests carried out on Rotavirus and Adenovirus-positive viral culture samples were positive,

with the development of three coloured lines. All 15 tests carried out on dilution buffer were negative, with the development of a single coloured line (control line).

To check inter-batch accuracy, one Rotavirus and Adenovirus-positive viral culture sample and a dilution buffer were tested 3 times on 3 different Combi-kit batches. All results were correct as expected. The 3 batches gave positive results for the Rotavirus and Adenovirus-positive viral culture sample and negative results for the dilution buffer for every test.

**D. Interference**

Cross-reactivity to samples positive for the following pathogens was tested and found to be negative: *Giardia lamblia*, *Escherichia coli* 0157: H7, *Salmonella typhimurium*, *Salmonella enteritidis*, *Escherichia coli* K99, *Coronavirus*, *Entamoeba histolytica*, *Entamoeba dispar*, *Cryptosporidium parvum*.

**XII. LIMITS OF THE KIT**

The test is qualitative and cannot predict the quantity of antigens present in the sample. Clinical presentation and other test results must be taken into consideration to establish diagnosis.

A positive test does not rule out the possibility that other pathogens may be present.

Kit test is an acute-phase screening test. Specimens that are collected after this phase may contain antigen titres below the reagent's sensitivity threshold. If a sample is given a negative result despite the observed symptoms, a culture should be started to check the sample.

**XIII. TECHNICAL PROBLEMS / COMPLAINTS**

If you encounter a technical problem or if performances do not correspond with those indicated in this package insert:

1. Record batch number of incriminated kit
2. If possible, keep problematic sample at cold for the time lapse of complaint management
3. Contact Coris BioConcept (info@corisbio.com) or your local distributor

**XIV. BIBLIOGRAPHIC REFERENCES**

- Th. Leclipteux, D. Col, M. Venuti, F. Paulart, D. Van Beers, M. De Foor and R. Viehof;** *Comparison of immunochromatography with ELISA to detect Adenovirus in stools specimens*; New Insights in gastrointestinal Diseases, London, UK, May, 1998
- R. Viehoff, D. Van Beers, M. De Foor, D. Col, M. Venuti, F. Paulart and Th. Leclipteux;** *Set up of a rapid immunochromatographic diagnostic test for Adenovirus detection*; European Society for Clinical Virology IV, Hamburg, Germany, August, 1998
- Sneyers et al.;** *Detection of rotavirus in faecal specimens with a monoclonal antibody enzyme-linked immunosorbent assay: comparison with polyclonal antibody enzyme-immunoassay and a latex agglutination test*; Comp. Immun. Microbiol. Infect. Dis., Vol 12, N° 4, pp 95-104, 1989
- E. Thomas et al.;** *Evaluation of Seven Immunoassays for Detection of Rotavirus in Pediatric Stool Samples*; J. Clin. Microbiol. 1988; 26 (6): 1189-1193
- Molyneaux et al.;** *Comparison of Six Commercial Kits for the Diagnosis of Rotavirus Infection in Man and Claves*; Serodiagn. and Immunother in Inf. Dis. 1989; 3: 123-134
- H. Dennehy, D.R. Gauntlett and W. E. Tente;** *Comparison of Nine Commercial Immunoassays for the Detection of Rotavirus in Fecal Specimens*; J. Clin. Microbiol. 1988; 26 (9): 1630-1634 - 6.E. & C. Kurstak
- Van Beers, M. De Foor, R. Viehoff, D. Col, M. Venuti and T. Leclipteux;** *Set-up of a new rapid immunochromatographic diagnostic test for Rotavirus detection*; Congress in Clinical Virology, Bologna, Settembre 1997, pag. 79

Last update: JANUARY 2011

<b>REF</b>	Catalogue number		Manufactured by
<b>IVD</b>	In vitro diagnostic medical device		Temperature limitation
	Contains sufficient for <n> tests	<b>DIL SPE</b>	Diluent specimen
	Consult instructions for use		Do not reuse
	Keep dry		Use by
<b>DIL AS</b>	Diluent assay	<b>CONT NaN<sub>3</sub></b>	Contains Sodium azide

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## XIII. PROBLEMES TECHNIQUES / RECLAMATIONS

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## XIV. REFERENCES BIBLIOGRAPHIQUES

1. **Th. Leclipteux, D. Col, M. Venuti, F. Paulart, D. Van Beers, M. De Foor and R. Viehoff;** Comparison of immunochromatography with ELISA to detect Adenovirus in stools specimens; New Insights in gastrointestinal Diseases, London, UK, May, 1998
2. **R. Viehoff, D. Van Beers, M. De Foor, D. Col, M. Venuti, F. Paulart and Th. Leclipteux;** Set up of a rapid immunochromatographic diagnostic test for Adenovirus detection; European Society for Clinical Virology IV, Hamburg, Germany, August, 1998
3. **Sneyers et al.;** Detection of rotavirus in faecal specimens with a monoclonal antibody enzyme-linked immunosorbent assay: comparison with polyclonal antibody enzyme-immunoassay and a latex agglutination test; Comp. Immun. Microbiol. Infect. Dis., Vol 12, N° 4, pp 95-104, 1989
4. **E. Thomas et al.;** Evaluation of Seven Immunoassays for Detection of Rotavirus in Pediatric Stool Samples; J. Clin. Microbiol. 1988; 26 (6): 1189-1193
5. **Molyneux et al.;** Comparison of Six Commercial Kits for the Diagnosis of Rotavirus Infection in Man and Claves; Serodiagn. and Immunother in Inf. Dis. 1989; 3: 123-134
6. **H. Dennehy, D.R. Gauntlett and W. E. Tente;** Comparison of Nine Commercial Immunoassays for the Detection of Rotavirus in Fecal Specimens; J. Clin. Microbiol. 1988; 26 (9): 1630-1634 - 6.E. & C. Kurstak
7. **Van Beers, M. De Foor, R. Viehoff, D. Col, M. Venuti and T. Leclipteux;** Set-up of a new rapid immunochromatographic diagnostic test for Rotavirus detection; Congress in Clinical Virology, Bologna, Settembre 1997, pag. 79

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