

# APT Agar

All purpose medium with Tween® proposed by EVANS and NIVEN (1951) and DEIBEL, EVANS and NIVEN (1957) for counting and cultivating heterofermentative lactic acid bacteria including *Lactobacillus*, *Leuconostoc* species, *Lactococcus lactis* and other microorganisms which require a high thiamine concentration in meat products, tinned foods, fruit juices and other foodstuffs.

The medium complies with the recommendations of the American Public Health Association (1992).

## Mode of Action

This medium contains a rich nutrient base with additives of Tween®, thiamine and several essential elements, which provide optimal growth conditions for the above mentioned bacteria. The culture medium is not selective, accompanying bacteria therefore also grow very well.

## Typical Composition (g/litre)

Peptone from casein 12.5; yeast extract 7.5; D(+)-glucose 10.0; sodium chloride 5.0; tri-sodium citrate 5.0; di-potassium hydrogen phosphate 5.0; Tween® 80 0.2; magnesium sulfate 0.8; manganese chloride 0.14; iron(II) sulfate 0.04; thiaminium dichloride 0.001; agar-agar 13.5.

## Preparation

Suspend 59.5 g/litre, fill into suitable containers, autoclave (15min at 121 °C). **Do not overheat**

pH: 6.7 ± 0.2 at 25 °C.

The prepared medium is clear and brown.

## Experimental Procedure and Evaluation

When performing bacterial counts, dilute the sample material and inoculate the APT Agar by the pour-plate method.

Incubation: 2 days at 35 °C aerobically.

In order to identify lactic acid bacteria which produce a green colouration, inoculate with the suspect colonies. After incubating for 24 hours at 32 °C, transfer a sample from the culture that has developed onto the cut surface of a smoked sausage. Place the sausage in a Petridish containing a damp piece of filter paper ("moist chamber"). Incubate for 18-24 hours at 32 °C and see whether there is a green colouration. A sample of the sausage which has not been inoculated serves as a control. In order to exclude other pigment-forming bacteria (e.g. *Pseudomonas*), a confirmatory bacteriological test (e.g. Gram-positive rods, negative catalase test, negative nitratase test, positive peroxidase test, acetoin production from glucose, ammonia production from arginine etc.) should also be performed.

## Literature

American Public Health Association: Compendium of Methods for the Microbiological Examination of Foods. - 3rd ed., 1992.

DEIBEL, R.H., EVANS, J.B., a. NIVEN, C.F.: Microbiological assay for the thiamin using *Lactobacillus viridescens*. - *J. Bact.*, **74**; 818-821 (1957).

EVANS J.B., a. NIVEN, C.F.: Nutrition of the heterofermentative *Lactobacilli* that cause greening of cured meat products. - *J. Bact.*, **62**; 599-603 (1951).

## Ordering Information

Product	Merck Cat. No.	Pack size
APT Agar	1.10453.0500	500 g

## Quality control

Test strains	Growth
<i>Lactobacillus acidophilus</i> ATCC 4356	good / very good
<i>Lactobacillus casei</i> ATCC 393	good / very good
<i>Lactobacillus fermentum</i> ATCC 9338	good / very good
<i>Lactobacillus plantarum</i> ATCC 14917	good / very good
<i>Lactobacillus viridescens</i> ATCC 12706	good / very good
<i>Leuconostoc mesenteroides</i> ATCC 9135	good / very good
<i>Lactococcus lactis</i> spp. <i>lactis</i> ATCC 19435	good / very good