



ChromoCult® Enterococci Agar

For detection of enterococci in water,
foodstuffs and other materials.



ChromoCult® Enterococci Agar



Ord. No. 1.00950.0500 (500 g)

Selective culture medium for the isolation, differentiation and enumeration of enterococci in water, foodstuffs and other materials. The medium contains specific chromogenic substrates – Enterococci grow as red colonies and are easily differentiated from other colony types.

Validated acc. to ISO 17994 : 2004 and approved by German EPA as alternativ method.

Mode of action

The presence of enterococci, especially *E. faecalis*, *E. faecium*, *E. durans* und *E. hirae* serves as an indicator for faecal contamination.

Growth of enterococci is stimulated by carefully selected peptones as well as phosphate and the addition of Tween® 80. Enterococci cleave the unique chromogenic substrates in the medium resulting in red colonies allowing a simple detection. Sodium azide and ox bile inhibit most accompanying microbial flora. Non-enterococci produce colourless, blue/violet or turquoise colonies. These colonies are easily distinguished from the red coloured colonies enterococci produce.

Typical composition (g/Litre)

Peptones 10.0; sodium chloride 5.0; sodium azide 0.2;
dipotassium hydrogenphosphate 3.4;
potassium dihydrogenphosphate 1.6; ox bile 0.5;
Tween® 80 1.0; chromogenic-mixture 0.25; Agar-agar 11.0

Preparation

Suspend 33.0 g in 1 litre of demin. water by heating in a boiling water bath or in a flowing stream of steam. Stir the contents to assist dissolution (approx. 45 minutes), let the medium cool to 45–50 °C. Pour into plates.

Do not autoclave! Do not overheat!

pH: 7.0 ±0.2 at 25 °C.

The plates are clear and slightly yellow. If stored at (+4 ±2 °C) and protected from light the plates are stabile for 2 weeks.

Experimental procedure

Inoculate the medium by the pour-plate method or by spreading the sample material on the surface of the plates. In addition the membrane-filter technique can also be used. The type of membrane filter affects the performance of the medium (growth and colouration of colonies). Best results are obtained using membrane filters of cellulose-mixed-ester material, e.g. Pall GN 6 (OSSMER, 1999)

Incubation: 24 ±4 hours at 35–37 °C.

Incubate for up to 44 ±4 hours to confirm a negative result.

Evaluation

Enterococci: red colonies with a diameter of 0.5 bis 2 mm.

Non-Enterococci: colourless (e.g. *Aerococcus viridans* ATCC 29503), blue/violet (e.g. *Aerococcus viridans* ATCC 10400), turquoise (e.g. *Streptococcus equi* ATCC 33398).





Colour makes the difference.


Quality control

Test strains	Inoculum (cfu/ml)	Growth	Colony colour
Enterococcus faecalis ATCC 19433	30 – 300	good	red
Enterococcus faecium ATCC 882	30 – 300	good	red
Enterococcus durans ATCC 6056	30 – 300	good	red
Enterococcus hirae ATCC 8043	30 – 300	good	red
Aerococcus viridans ATCC 10400	1000 – 2000	fair/none	blue/violet
Bacillus cereus ATCC 11778	1000 – 2000	-	-
Escherichia coli ATCC 11775	-	-	-
Pseudomonas aeruginosa ATCC 27853	-	-	-

Aerococcus viridans



Enterococci



Literature

ALTHAUS, H.; DOTT, W.; HAVEMEISTER, G.; MÜLLER, H. E. and SACRÉ, C.:
Faecal streptococci as indicator organisms of drinking water.

• *Zbl. Bakt. Hyg., I. Abt. Orig. A* 252: 154–165 (1982)

OSSMER, R.; SCHMIDT, W.; MENDE, U.:

ChromoCult® Coliform Agar – Influence of Membrane Filter Quality on Performance. Poster presentation Congresso de la Sociedad. Espanola de Microbiologia, Granada, Spain

• (1999)

AMOROS, I.:

Evaluation of ChromoCult® Enterococci Broth (with Agar).

Poster presentation Congress of Spanish Society of Microbiology. Madrid

• (1995)

LITSKY, W.; MALLMANN, W. L. and FIFIELD, C. W.:

A new medium for the detection of enterococci in water.

• *Amer. J. Pbl. Hlth.* 43: 873–879 (1953)

MANAFI, M. and WINDHAGER, K.:

Rapid identification of enterococci in water with a new chromogenic assay.

• *Abstr. P-107, pp. 453, Abstracts of the 97th Meeting of the American Society for Microbiology, Miami, USA (1997)*

SNYDER, M. L. and LICHSTEIN, H. C.:

Sodium azide as an inhibiting substance for Gram-negative bacteria.

• *J. Infect. Dis.* 67: 113–115 (1940)



W 286111 11/06

For further information on Merck
and our products contact

Merck KGaA
64271 Darmstadt, Germany
Fax: 0049 (0)61 51 72 - 60 80
E-mail: microbiology@merck.de
www.merck.de
microbiology.merck.de

We provide information and advice to our customers to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose. ChromoCult® is a registered trademark of Merck KGaA, Darmstadt, Germany.