**remel**

# ***RapID™ Yeast Plus System***

**English - p 2**

For the biochemical identification of medically important yeast and yeast-like organisms

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
## RapID™ Yeast Plus System (English)

### INTENDED USE

Remel RapID™ Yeast Plus System is a qualitative micromethod employing conventional and chromogenic substrates for the identification of medically important yeast, yeast-like, and related organisms isolated from human clinical specimens. A complete listing of the organisms addressed by the RapID™ Yeast Plus System is given in the RapID™ Yeast Plus Differential Chart.

### SUMMARY AND EXPLANATION

The RapID™ Yeast Plus System is comprised of (1) RapID™ Yeast Plus Panels, (2) RapID™ Yeast Plus Reagent A and (3) RapID™ Yeast Plus Reagent B. Each RapID™ Yeast Plus Panel has several reaction cavities molded into the periphery of a plastic disposable tray. Reaction cavities contain dehydrated reactants and the tray allows the simultaneous inoculation of each cavity with a predetermined amount of inoculum. A suspension of the test organism in RapID™ Inoculation Fluid is used as the inoculum which rehydrates and initiates test reactions. After incubation of the panel, each test cavity is examined for reactivity by noting the development of a color. In some cases, reagents must be added to the test cavities to provide a color change. The resulting pattern of positive and negative test scores is used as the basis for identification of the test isolate by comparison of test results to reactivity patterns stored in a database or through the use of a computer-generated Code Compendium.

### PRINCIPLE

The tests used in the RapID™ Yeast Plus System are based upon the microbial degradation of specific substrates detected by various indicator systems. The reactions employed are a combination of conventional tests and single-substrate chromogenic tests and are described below in Table 1.

### REAGENTS\*

#### RapID™ Yeast Plus Reagent A (provided with kit) (15 ml/Btl.)

Reactive Ingredient per Liter:  
Potassium Hydroxide.....16.0 g

#### RapID™ Yeast Plus Reagent B (provided with kit) (10 ml/Btl.)

Reactive Ingredient per Liter:  
p-Dimethylaminocinnamaldehyde.....0.06 g

#### RapID™ Inoculation Fluid (REF 8325106 supplied separately) (2 ml/Tube)

KCl.....6.0 g  
CaCl<sub>2</sub>.....0.5 g  
Demineralized Water.....1000.0 ml

\*Adjusted as required to meet performance standards.

### PRECAUTIONS

This product is For *In Vitro* Diagnostic Use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, media, and test panels after their use. Directions should be read and followed carefully.

**Table 1. Principles and Components of the RapID™ Yeast Plus System**

Cavity #	Test Code	Reactive Ingredient	Quantity	Principle	Bibliography #
1	GLU	Glucose	1.0%	Utilization of the carbohydrate produces acidic products which lower the pH and change the indicator.	1
2	MAL	Maltose	1.0%		
3	SUC	Sucrose	1.0%		
4	TRE	Trehalose	1.0%		
5	RAF	Raffinose	1.0%		
6	LIP	Fatty acid ester	1.0%	Hydrolysis of the fatty acid ester releases acidic products which lower the pH and change the indicator.	2
7	NAGA	p-Nitrophenyl-N-acetyl-β,D-galactosaminide	0.05%	Enzymatic hydrolysis of the colorless aryl-substituted glycoside or phosphoester releases yellow o- or p-nitrophenyl which is detected with RapID™ Yeast Plus Reagent A.	3-8
8	αGLU	p-Nitrophenyl-α,D-glucoside	0.05%		
9	βGLU	p-Nitrophenyl-β,D-glucoside	0.05%		
10	ONPG	o-Nitrophenyl-β,D-galactoside	0.05%		
11	αGAL	p-Nitrophenyl-α,D-galactoside	0.05%		
12	βFUC	p-Nitrophenyl-β,D-fucoside	0.05%		
13	PHS	p-Nitrophenyl phosphate	0.05%		
14	PCHO	p-Nitrophenyl phosphorylcholine	0.05%		
15	URE	Urea	0.3%	Hydrolysis of urea produces basic products which raise the pH and change the indicator.	9
16	PRO	Proline-β-naphthylamide	0.01%	Enzymatic hydrolysis of the arylamide substrate releases free β-naphthylamine which is detected with RapID™ Yeast Plus Reagent B	10
17	HIST	Histidine β-naphthylamide	0.01%		
18	LGY	Leucyl-glycine β-naphthylamide	0.01%		

### Caution!

1. RapID™ Yeast Plus Reagent A may cause irritation to skin, eyes, and respiratory system.
2. RapID™ Yeast Plus Reagent B is toxic and may cause harm to the environment. Harmful by inhalation, contact with skin or eyes, or if swallowed. May impair fertility or cause harm to unborn child.
3. Refer to Material Safety Data Sheet for detailed information on reagent chemicals.

### STORAGE

The RapID™ Yeast Plus System should be stored in its original container at 2-8°C until used. Allow product to come to room temperature before use. DO NOT interchange reagents among different RapID™ systems. Remove only the number of panels necessary for testing. Reseal the plastic pouch and promptly return to 2-8°C. Panels must be used the same day they are removed from storage. RapID™ Inoculation Fluid should be stored in its original container at room temperature (20-25°C) until used.

### PRODUCT DETERIORATION

This product should not be used if (1) the color of the reagent has changed, (2) the expiration date has passed, (3) the plastic tray is broken or the lid is compromised, or (4) there are other signs of deterioration.

### SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.<sup>11,12</sup>

### MATERIALS SUPPLIED

(1) 20 RapID™ Yeast Plus Panels, (2) 20 report forms, (3) 1 each, RapID™ Yeast Plus Reagents A and B (plastic dropper-bottles containing reagent sufficient for 20 panels), (4) 2 chipboard incubation trays, (5) 1 RapID™ Yeast Plus Inoculation Card, (6) Instructions for use.

### MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swab, collection containers, (3) Incubators, alternative environmental systems, (4) Supplemental media, (5) Quality control organisms, (6) Gram stain reagents or demineralized water, (7) Microscope slides, (8) Cotton swabs, (9) RapID™ Inoculation Fluid - 2 ml (REF 8325106), (10) Pipettes, (11) RapID™ Yeast Plus Code Compendium (REF 8327007) or ERIC® (Electronic RapID™ Compendium, REF 8323600).

### PROCEDURE

#### Inoculum Preparation

1. Test organisms must be grown in pure culture and examined by Gram stain or wet mount prior to use in the system.

**Note:** Only organisms which demonstrate yeast-like appearance and growth characteristics should be tested using the RapID™ Yeast Plus System.

2. The following medium is recommended:  
Sabouraud Dextrose Agar (SDA) - Emmons formulation

#### Notes:

- Plates used for inoculum preparation should be incubated at 30°C and should be 48 hours old.
- The use of media other than that recommended may compromise test performance.


3. Using a cotton swab or inoculating loop, suspend sufficient growth from the agar plate culture in RapID™ Inoculation Fluid (2 ml) to achieve a visual turbidity as detailed under **Notes**, using the RapID™ Yeast Plus Inoculation Card.

### Notes:


- Select well-isolated colonies of the test isolate and add **incrementally** to the RapID™ Inoculation Fluid to avoid clumping and over-inoculation. Continue to add organism until the turbidity of the suspension **completely** obliterates the black lines on the Inoculation Card. Once the black lines on the Inoculation Card are no longer visible, inoculum preparation is complete.
  - Suspensions significantly less turbid than the required inoculum density will result in aberrant reactions.
  - Suspensions that are **slightly** more turbid than the required inoculum density will not affect test performance and are recommended for stock cultures and quality control strains. However, suspensions that are **significantly** more turbid will compromise test performance.
  - Suspensions should be mixed thoroughly and vortexed if required.
  - Suspensions should be used within 15 minutes of preparation.
3. An agar plate may be inoculated for purity and any additional testing that may be required using a loopful of the test suspension from the inoculation fluid tube. Incubate the plate for 24-72 hours at 30°C.

### Inoculation of RapID™ Yeast Plus Panels

1. Peel back the lid of the panel over the inoculation port by pulling the tab marked "Peel to Inoculate" up and to the left.
2. Using a pipette, gently transfer the **entire** contents of the Inoculation Fluid tube into the upper right-hand corner of the panel. Reseal the inoculation port of the panel by pressing the peel-back tab back in place.
3. After adding the test suspension, and while keeping the panel on a level surface, tilt the panel back away from the reaction cavities at approximately a 45-degree angle (see below).


4. While tilted back, gently rock the panel from side to side to evenly distribute the inoculum along the rear baffles as illustrated below.

5. While maintaining a level, horizontal position (best achieved by using the bench top against the reaction cavity bottoms), slowly tilt the panel forward toward the reaction cavities until the inoculum flows along the baffles into the reaction cavities (see below). This should evacuate all of the inoculum from the rear portion of the panel.

**Note:** If the panel is tilted too quickly, air may be trapped at the test cavity junction, restricting fluid movement.

- Return the panel to a level position. If necessary, gently tap the panel on the bench top to remove any air trapped in the cavities.

**Notes:**

- Examine the test cavities, which should appear bubble-free and uniformly filled. Slight irregularities in test cavity fills are acceptable and will not affect test performance. If the panel is grossly misfilled, a new panel should be inoculated and the misfilled panel discarded.
- Complete the inoculation of each panel receiving inoculation fluid before inoculating additional panels.
- Do not allow the inoculum to rest in the back portion of the panel for prolonged periods without completing the procedure.

**Incubation of RapID™ Yeast Plus Panels**

Incubate inoculated panels at 30°C in a non-CO<sub>2</sub> incubator for 4 hours. For ease of handling, panels may be incubated in the chipboard incubation trays provided with the kit.

**Notes:**

- If a 30°C, non-CO<sub>2</sub> incubator is unavailable, RapID™ Yeast Plus Panels may be incubated at controlled room temperature.
- Incubation of RapID™ Yeast Plus Panels at 35-37°C may result in aberrant reactions.

**Scoring of RapID™ Yeast Plus Panels**

RapID™ Yeast Plus panels contain 18 reaction cavities that provide 18 test scores. Reagent-requiring tests (cavities 7-14 and 16-18) are designated with a box drawn around them.

- While firmly holding the RapID™ Yeast Plus panel on the benchtop, peel off the label lid over the reaction cavities by pulling the lower right hand tab up and to the left.
- Add the following reagents to the cavities indicated:
  - Add 1 drop of RapID™ Yeast Plus Reagent A to cavities 7 (NAGA) through 14 (PCHO).
  - Add 1 drop of RapID™ Yeast Plus Reagent B to cavities 16 (PRO) through 18 (LGY).
- After the addition of RapID™ Yeast Plus Reagent B, allow at least 30 seconds but no more than 1 minute for color development.
 

**Note:** Cavities exhibiting layers of color may be mixed using an applicator stick prior to reading.
- Read and score the test cavities from left to right using the interpretation guide presented in Table 2. Record the scores in the appropriate boxes on the report form.
- Reference the microcode obtained on the report form in the RapID™ Yeast Plus Code Compendium or ERIC® (Electronic RapID™ Compendium) for the identification.

**RapID™ Yeast Plus Panel Test Location**

Cavity #	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
Test Code	GLU	MAL	SUC	TRE	RAF	LIP	NAGA	αGLU	βGLU	ONPG	αGAL	βFUC	PHS	PCHO	URE	PRO	HIST	LGY
							RapID™ Yeast Plus Reagent A								RapID™ Yeast Plus Reagent B			

**Table 2. Interpretation of RapID™ Yeast Plus System Tests\***

Cavity #	Test Code	Reagent	Reaction		Comments
			Positive	Negative	
1	GLU				
2	MAL				
3	SUC	None	Yellow	Blue, green-blue, or green	Only the development of a distinct yellow color should be scored as positive.
4	TRE				
5	RAF				
6	LIP	None	Yellow	Red, pink, orange, or gold	Only the development of a distinct lemon-yellow color should be scored as positive.
7	NAGA				
8	αGLU				
9	βGLU				
10	ONPG	RapID™ Yeast Plus Reagent A	Yellow	Clear or cream button	Any shade of yellow color development should be scored as positive.
11	αGAL				
12	βFUC				
13	PHS				
14	PCHO				
15	URE	None	Red or dark red-orange	Yellow, yellow-orange, or orange	Only the development of a red or dark red-orange color should be scored as positive. Other shades of orange should be scored as negative.
16	PRO				
17	HIST	RapID™ Yeast Plus Reagent B	Purple, red, or dark pink	Clear, straw, orange, or pale to medium pink	Only the development of a distinct purple, red, or dark pink color should be scored as positive. Pale shades should be scored as negative.
18	LGY				

\*NOTE: Panels should be read by looking down through the reaction cavities against a white background.

## RESULTS AND RANGE OF EXPECTED VALUES

The RapID™ Yeast Plus Differential Chart illustrates the expected results for the RapID™ Yeast Plus System. Differential Chart results are expressed as a series of positive percentages for each system test. This information statistically supports the use of each test and provides the basis, through numerical coding of digital test results, for a probabilistic approach to the identification of the test isolate.

Identifications are made using individual test scores from RapID™ Yeast Plus panels in conjunction with other laboratory information to produce a pattern that statistically resembles known reactivity for taxa recorded in the RapID™ System database. These patterns are compared through the use of the RapID™ Yeast Plus Differential Chart, or by derivation of a microcode and the use of the RapID™ Yeast Plus Code Compendium or ERIC®.

## QUALITY CONTROL

All lot numbers of the RapID™ Yeast Plus System have been tested using the following quality control organisms and have been found to be acceptable. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not

be reported. Table 3 lists expected results for the selected battery of test organisms.

### Notes:

- RapID™ Yeast Plus Reagent quality control is accomplished by obtaining the expected reactions for tests requiring the addition of the reagents (cavities 7-14 for Reagent A; cavities 16-18 for Reagent B).
- Organisms which have been repeatedly transferred on agar media for prolonged periods may provide aberrant results.
- Quality control strains should be stored frozen or lyophilized, or on SDA slants (Emmons formulation) at 2-8°C. Prior to use, quality control strains should be transferred 2-3 times from storage to SDA (Emmons formulation). The final subculture to be used for QC testing should be incubated at 30°C for 48 hours.
- Formulations, additives, and ingredients of culture media vary from manufacturer to manufacturer and may vary from batch to batch. As a result, culture media may influence constitutive enzymatic activity of designated quality control strains. If quality control strain results differ from the patterns indicated, a subculture onto medium from a different batch or from another manufacturer will often resolve quality control discrepancies.

**Table 3. Quality Control Chart for RapID™ Yeast Plus Panels**

Organism	GLU	MAL	SUC	TRE	RAF	LIP	NAGA	αGLU	βGLU	ONPG	αGAL	βFUC	PHS	PCHO	URE	PRO	HIST	LGY
<i>Candida albicans</i> ATCC® 14053	+	+	-	-	-	-	+	+	-	-	-	-	V	-	-	+	V	V
<i>Candida glabrata</i> ATCC® 2001	+	-	-	+	-	V	-	-	-	-	-	-	-	-	-	-	V	V
<i>Candida kefyr</i> ATCC® 2512	+	V	+	V	+	V	-	-	+	+	-	+	-	-	-	V	V	V
<i>Cryptococcus laurentii</i> ATCC® 66036	-	-	-	-	-	V	-	+	+	-	+	V	V	V	+	-	-	-
<i>Yarrowia lipolytica</i> ATCC® 9773	-	-	-	-	-	+	-	-	V	-	-	-	+	+	-	+	+	+

+, positive; -, negative; V, variable

## LIMITATIONS

1. The use of the RapID™ Yeast Plus System and the interpretation of results require a competent laboratorian who is trained in general microbiological methods and who should judiciously make use of knowledge, experience, specimen information, and other pertinent procedures before reporting the isolate identity obtained using the RapID™ Yeast Plus System.
2. The RapID™ Yeast Plus System must be used with pure cultures of test organisms. The use of mixed microbial populations or direct testing of clinical material without culture will result in aberrant results.
3. The RapID™ Yeast Plus System is designed for use with the taxa listed in the RapID™ Yeast Plus Differential Chart. The use of organisms not specifically listed may lead to misidentifications.
4. Expected values listed for RapID™ Yeast Plus System tests may differ from conventional test results or previously reported information.
5. The accuracy of the RapID™ Yeast Plus System is based upon the statistical use of a multiplicity of specially designed tests and an exclusive, proprietary database. The use of any single

test found in the RapID™ Yeast Plus System to establish the identification of a test isolate is subject to the error inherent in that test alone.

## PERFORMANCE CHARACTERISTICS

The RapID™ Yeast Plus System performance characteristics have been established by laboratory testing of 500 clinical, reference, and type cultures at Remel. Overall, the RapID™ Yeast Plus System correctly identified 476 (95.2%) of the organisms tested. A total of 378 isolates were compared using the RapID™ Yeast Plus and the API 20C.<sup>13</sup> The RapID™ Yeast Plus System agreed with the API 20C for 361 (95.5%) of the isolates tested.

The RapID™ Yeast Plus System has been independently evaluated using 185 clinical yeast isolates.<sup>14</sup> A total of 181 (97.8%) isolates were correctly identified by the RapID™ Yeast Plus System without additional testing, and 4 isolates (2.2%) were correctly identified after additional testing was performed. No misidentifications were observed.

## RAPID™ Yeast Plus Differential Chart

Organism	GLU	MAL	SUC	TRE	RAF	LIP	NAGA	αGLU	βGLU	ONPG	αGAL	βFUC	PHS	PCHO	URE	PRO	HIST	LGY
<i>Aureobasidium pullulans</i>	29	4	19	5	0	74	84	81	80	1	81	7	92	78	96	95	18	90
<i>Blastoschizomyces capitatus</i>	16	9	8	6	4	74	0	21	41	0	0	0	5	2	0	5	98	96
<i>Candida albicans</i>	96	94	2	14	1	3	90	94	4	0	0	0	76	2	1	99	36	9
<i>C. apicola</i>	98	0	94	1	1	0	0	1	1	0	2	7	2	9	0	0	96	95
<i>C. ciferrii</i>	3	0	2	0	0	58	96	79	94	0	11	2	0	86	2	0	55	11
<i>C. colliculosa</i>	98	5	98	97	1	0	0	95	90	1	0	2	89	0	3	99	86	66
<i>C. famata</i> <sup>1</sup>	90	0	77	4	8	0	0	49	74	7	7	2	76	26	0	99	18	4
<i>C. glabrata</i> <sup>2</sup>	98	8	2	96	2	24	1	4	11	0	1	0	1	1	0	0	29	13
<i>C. guilliermondii</i>	91	0	92	4	38	3	0	63	40	0	70	1	91	88	7	99	70	29
<i>C. intermedia</i>	98	44	96	74	72	1	0	33	99	0	0	0	96	32	0	96	97	28
<i>C. kefyi</i> <sup>3</sup>	97	7	98	3	92	6	0	4	93	73	0	88	0	0	3	1	16	4
<i>C. krusei</i>	99	1	1	1	0	1	0	0	26	0	0	0	1	0	9	0	69	35
<i>C. lambica</i>	96	0	0	0	0	4	0	0	8	0	0	0	90	91	1	92	78	8
<i>C. lusitanae</i>	99	0	10	8	0	2	0	97	98	0	0	0	90	88	6	99	66	6
<i>C. marina</i>	95	0	0	0	0	0	0	1	98	95	0	92	88	90	2	98	92	2
<i>C. parapsilosis</i>	98	3	4	0	0	2	0	94	9	0	0	1	80	77	4	98	15	6
<i>C. rugosa</i>	0	0	0	0	0	84	14	0	2	60	0	0	16	5	2	16	77	9
<i>C. stellatoidea</i>	98	95	1	9	2	9	95	95	6	1	1	1	81	2	2	0	42	9
<i>C. tropicalis</i> <sup>4</sup>	98	64	87	84	3	4	14	95	5	0	1	1	60	67	2	1	31	28
<i>C. utilis</i>	98	1	96	1	96	4	0	2	90	0	1	5	88	92	2	17	90	85
<i>C. zeylanoides</i>	7	0	0	0	0	0	0	0	20	0	0	0	12	2	1	98	2	1
<i>Cryptococcus albidus</i>	86	61	81	69	11	14	0	90	92	0	1	18	90	88	90	70	61	11
<i>Cr. humicolus</i> <sup>5</sup>	8	3	0	0	2	1	38	86	90	0	98	66	88	93	81	65	17	36
<i>Cr. laurentii</i>	4	0	1	0	0	0	1	81	77	0	91	16	88	85	98	8	5	4
<i>Cr. neoformans</i>	68	16	44	5	11	3	15	12	36	0	0	0	11	1	98	1	2	2
<i>Cr. terreus</i>	0	0	0	0	0	0	0	0	19	0	98	0	74	70	70	98	66	21
<i>Cr. uniguttulatus</i>	2	2	2	2	0	2	1	1	96	2	8	13	23	19	99	95	25	39
<i>Geotrichum</i> spp.	2	0	0	0	0	2	0	0	84	0	0	71	1	0	1	0	90	61
<i>Hanseniaspora guilliermondii</i> / <i>uvarum</i>	98	0	0	0	5	2	0	2	98	0	0	81	0	1	0	1	16	2
<i>Hansenula wingei</i>	98	0	0	0	0	2	1	90	97	0	0	90	90	93	0	13	90	24
<i>Kluyveromyces</i> spp.	98	2	98	2	96	2	0	13	7	0	0	0	6	2	0	55	3	37
<i>Pichia anomala</i> <sup>6</sup>	99	22	96	1	94	2	0	90	98	0	0	82	86	88	2	1	31	7
<i>Prototheca wickerhamii</i>	98	76	85	92	74	64	0	1	0	0	0	0	66	90	1	90	8	3
<i>P. zopfii</i>	98	70	80	63	71	8	0	0	2	0	0	1	90	1	0	3	15	5
<i>Rhodotorula glutinis</i>	1	1	1	1	1	33	0	34	0	0	1	33	2	1	99	96	5	3
<i>R. minuta</i>	2	2	2	2	1	38	3	90	81	0	2	88	23	2	98	29	7	21
<i>R. rubra</i>	9	6	5	3	2	46	0	64	2	0	0	0	28	26	97	99	71	62
<i>Saccharomyces cerevisiae</i>	98	24	96	1	92	2	1	85	82	5	9	1	3	1	1	1	66	44
<i>Sporobolomyces salmonicolor</i>	0	0	0	0	0	0	11	0	0	0	0	0	0	0	97	98	88	96
<i>Trichosporon beigelii</i>	5	0	1	1	0	47	79	46	93	0	71	58	60	64	78	43	47	62
<i>Yarrowia lipolytica</i>	0	0	0	0	0	95	0	0	32	0	0	1	73	61	3	98	90	86

<sup>1</sup>Previously designated as *Torulopsis candida*

<sup>2</sup>Previously designated as *Torulopsis glabrata*

<sup>3</sup>Previously designated as *Candida pseudotropicalis*

<sup>4</sup>Includes species previously designated as *Candida paratropicalis*

<sup>5</sup>Previously designated as *Candida humicola*

<sup>6</sup>Previously designated as *Hansenula anomala*

**BIBLIOGRAPHY**

1. Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy. 1991. Manual of Clinical Microbiology. 5<sup>th</sup> ed., p. 617-629. ASM, Washington, D.C.
2. Lodder, J. 1970. The Yeasts - A Taxonomic Study. North Holland Publishing Co., Amsterdam - London.
3. Bobey, D.G., J.J. Bradna, and D.B. Florek-Ebeling. 1980. Abstract C-254. Abstracts of the 80<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
4. Bobey, D.G. and G.M. Ederer. 1981. J. Clin. Microbiol. 13:393-394.
5. David, H.L. and M.T. Jahan. 1977. J. Clin. Microbiol. 5:383-384.
6. Perry, J.L., G.R. Miller, and D.L. Carr. 1990. J. Clin. Microbiol. 28:614-615.
7. Smith, R.F., D. Blasi, and S.L. Dayton. 1973. Appl. Microbiol. 26:364-367.
8. Smitka, C.M. and S.G. Jackson. 1989. J. Clin. Microbiol. 27:203-206.
9. Roberts, G.D., C.D. Horstmeier, G.A. Land, and J.H. Foxworth. 1978. J. Clin. Microbiol. 7:584-588.
10. Norris, J.R. and D.W. Ribbons. 1976. Methods in Microbiology. Vol. 9, p. 1-14. Academic Press, New York, NY.
11. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenover. 2003. Manual of Clinical Microbiology. 8<sup>th</sup> ed., Vol. 2. ASM, Washington, D.C.
12. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. Bailey and Scott's Diagnostic Microbiology. 11<sup>th</sup> ed. Mosby, St. Louis, MO.
13. Marler, J.K. and L.A. Eriquez. 1995. Abstract C-418. Abstracts of the 95<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
14. Kitch, T. and P.C. Appelbaum. 1995. Abstract C-419. Abstracts of the 95<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
15. Eriquez, L.A. and J.K. Marler. 1994. Abstract C-490. Abstracts of the 94<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.

16. Lee, K.L., M.E. Reza, R.R. Watson, and C.C. Campbell. 1975. Sabouraudia 13:132-141.
17. Marler, J.K. and L.A. Eriquez. 1994. Abstract C-489. Abstracts of the 94<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.

**PACKAGING**

REF 8311007, RapID™ Yeast Plus System ..... 20 Tests/Kit

**Symbol Legend**


<b>REF</b>	Catalog Number
<b>IVD</b>	In Vitro Diagnostic Medical Device
<b>LAB</b>	For Laboratory Use
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
<b>LOT</b>	Batch Code (Lot Number)
	Use By (Expiration Date)
<b>EC REP</b>	European Authorized Representative

RapID™ is a trademark of Remel Inc.  
 ERIC® is a registered trademark of Remel Inc.  
 ATCC® is a registered trademark of American Type Culture Collection.

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1. Balows, A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy. 1991. Manual of Clinical Microbiology. 5<sup>th</sup> ed., p. 617-629. ASM, Washington, D.C.
2. Lodder, J. 1970. The Yeasts - A Taxonomic Study. North Holland Publishing Co., Amsterdam - London.
3. Bobey, D.G., J.J. Bradna, and D.B. Florek-Ebeling. 1980. Abstract C-254. Abstracts of the 80<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
4. Bobey, D.G. and G.M. Ederer. 1981. J. Clin. Microbiol. 13:393-394.
5. David, H.L. and M.T. Jahan. 1977. J. Clin. Microbiol. 5:383-384.
6. Perry, J.L., G.R. Miller, and D.L. Carr. 1990. J. Clin. Microbiol. 28:614-615.
7. Smith, R.F., D. Blasi, and S.L. Dayton. 1973. Appl. Microbiol. 26:364-367.
8. Smitka, C.M. and S.G. Jackson. 1989. J. Clin. Microbiol. 27:203-206.
9. Roberts, G.D., C.D. Horstmeier, G.A. Land, and J.H. Foxworth. 1978. J. Clin. Microbiol. 7:584-588.
10. Norris, J.R. and D.W. Ribbons. 1976. Methods in Microbiology. Vol. 9, p. 1-14. Academic Press, New York, NY.
11. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenover. 2003. Manual of Clinical Microbiology. 8<sup>th</sup> ed., Vol. 2. ASM, Washington, D.C.
12. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. Bailey and Scott's Diagnostic Microbiology. 11<sup>th</sup> ed. Mosby, St. Louis, MO.
13. Marler, J.K. and L.A. Eriquez. 1995. Abstract C-418. Abstracts of the 95<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
14. Kitch, T. and P.C. Appelbaum. 1995. Abstract C-419. Abstracts of the 95<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
15. Eriquez, L.A. and J.K. Marler. 1994. Abstract C-490. Abstracts of the 94<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
16. Lee, K.L., M.E. Reza, R.R. Watson, and C.C. Campbell. 1975. Sabouraudia 13:132-141.

17. Marler, J.K. and L.A. Eriquez. 1994. Abstract C-489. Abstracts of the 94<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.

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<sup>6</sup>Precedentemente denominato *Hansenula anomala*.

**BIBLIOGRAFIA**

1. Balows, A., W.J. Hausler, K.L. Herrmann, H.D. Isenberg, and H.J. Shadomy. 1991. Manual of Clinical Microbiology. 5<sup>th</sup> ed., p. 617-629. ASM, Washington, D.C.
2. Lodder, J. 1970. The Yeasts - A Taxonomic Study. North Holland Publishing Co., Amsterdam - London.
3. Bobey, D.G., J.J. Bradna, and D.B. Florek-Ebeling. 1980. Abstract C-254. Abstracts of the 80<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
4. Bobey, D.G. and G.M. Ederer. 1981. J. Clin. Microbiol. 13:393-394.
5. David, H.L. and M.T. Jahan. 1977. J. Clin. Microbiol. 5:383-384.
6. Perry, J.L., G.R. Miller, and D.L. Carr. 1990. J. Clin. Microbiol. 28:614-615.
7. Smith, R.F., D. Blasi, and S.L. Dayton. 1973. Appl. Microbiol. 26:364-367.
8. Smitka, C.M. and S.G. Jackson. 1989. J. Clin. Microbiol. 27:203-206.
9. Roberts, G.D., C.D. Horstmeier, G.A. Land, and J.H. Foxworth. 1978. J. Clin. Microbiol. 7:584-588.
10. Norris, J.R. and D.W. Ribbons. 1976. Methods in Microbiology. Vol. 9, p. 1-14. Academic Press, New York, NY.
11. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 2003. Manual of Clinical Microbiology. 8<sup>th</sup> ed., Vol. 2. ASM, Washington, D.C.
12. Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. Bailey and Scott's Diagnostic Microbiology. 11<sup>th</sup> ed. Mosby, St. Louis, MO.
13. Marler, J.K. and L.A. Eriquez. 1995. Abstract C-418. Abstracts of the 95<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
14. Kitch, T. and P.C. Appelbaum. 1995. Abstract C-419. Abstracts of the 95<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.
15. Eriquez, L.A. and J.K. Marler. 1994. Abstract C-490. Abstracts of the 94<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.

16. Lee, K.L., M.E. Reza, R.R. Watson, and C.C. Campbell. 1975. Sabouraudia 13:132-141.
17. Marler, J.K. and L.A. Eriquez. 1994. Abstract C-489. Abstracts of the 94<sup>th</sup> General Meeting of the American Society for Microbiology. ASM, Washington, D.C.

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RapID™ es una marca comercial de Remel Inc.
ERIC® es una marca registrada de Remel Inc.