remel

RapIDTM ***STAPH PLUS System***

English - p 2

For the biochemical identification of Staphylococci
and related organisms

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
RapID™ STAPH PLUS System (English)

INTENDED USE

Remel RapID™ STAPH PLUS System is a qualitative micromethod employing conventional and chromogenic substrates for the identification of medically important *Staphylococcus* species and related organisms that have been isolated from human clinical specimens. The RapID™ STAPH PLUS System is designed for use with microorganisms that routinely grow on plated media used for the isolation and cultivation of staphylococci and other gram-positive cocci in the clinical laboratory. A complete listing of the organisms addressed by the RapID™ STAPH PLUS System is given in the RapID™ STAPH PLUS Differential Chart.

SUMMARY AND EXPLANATION

The RapID™ STAPH PLUS System is comprised of (1) RapID™ STAPH PLUS Panels and (2) RapID™ STAPH PLUS Reagent. Each RapID™ STAPH PLUS Panel has 18 reaction cavities molded into the periphery of a plastic disposable tray. Reaction cavities contain dehydrated reactants and the tray allows the simultaneous inoculation of each cavity with a predetermined amount of inoculum. A suspension of the test organism in RapID™ Inoculation Fluid is used as the test inoculum, which rehydrates and initiates the test reactions. After incubation of the panel, each test cavity is examined for reactivity by noting the development of a color. In some cases, reagents must be added to the test cavities to provide a color change. The resulting pattern of positive and negative test scores is used as the basis for identification of the test isolate by comparison of results to reactivity patterns stored in a database through the use of the Electronic RapID™ Compendium (ERIC®) or by use of the RapID™ STAPH PLUS Differential Chart.

PRINCIPLE

The tests used in RapID™ STAPH PLUS System are based upon the microbial degradation of specific substrates detected by various indicator systems. The reactions employed are a combination of conventional and single-substrate chromogenic tests and are described in Table 1.

REAGENTS*

RapID™ STAPH PLUS Reagent (provided with kit) (15 ml/Btl)
 Reactive ingredient per liter:
 p-Dimethylaminocinnamaldehyde0.06 g

RapID™ STAPH PLUS Reagent is intended for use only in conjunction with RapID™ STAPH PLUS System panels.

RapID™ Inoculation Fluid (REF R8325106 supplied separately)
 (2 ml/Tube)

KCl 6.0 g
 CaCl₂ 0.5 g
 Demineralized Water 1000.0 ml

RapID™ Nitrate A Reagent (REF R8309003 supplied separately)
 (15 ml/Btl)

Sulfanilic Acid 8.0 g
 Glacial Acetic Acid 280.0 ml
 Demineralized Water 720.0 ml

RapID™ Nitrate B Reagent (REF R8309004 supplied separately)
 (15 ml/Btl)

n,n-Dimethyl-1-naphthylamine 6.0 g
 Glacial Acetic Acid 280.0 ml
 Demineralized Water 720.0 ml

*Adjusted as required to meet performance standards.

PRECAUTIONS

This product is For *In Vitro* Diagnostic Use and should be used by properly trained individuals. Precautions should be taken against the dangers of microbiological hazards by properly sterilizing specimens, containers, media, and test panels after use. Directions should be read and followed carefully.

Caution!

1. RapID™ Nitrate A Reagent and RapID™ Nitrate B Reagent may cause irritation to skin, eyes, and respiratory system.
2. RapID™ STAPH PLUS Reagent is toxic and may cause harm to the environment. Harmful by inhalation, contact with skin or eyes, or if swallowed. May impair fertility or cause harm to unborn child.
3. Refer to Material Safety Data Sheet for detailed information on reagent chemicals.

Table 1. Principles and Components of the RapID™ STAPH PLUS System

Cavity #	Test Code	Reactive Ingredient	Quantity	Principle	Bibliography #
1	ADH	Arginine	2%	Utilization of the amino acid substrate produces basic products, which raise the pH and change the indicator.	1,9,11,14,16
2	ODC	Ornithine	2%		
3	LIP	Fatty acid ester	2%	Utilization of the fatty acid substrate produces acidic products, which lower the pH and change the indicator.	1,9,16
4	SUC	Sucrose	2%	Utilization of the carbohydrate substrate produces acidic products, which lower the pH and change the indicator.	1,9,11,14,16
5	MANO	Mannose	2%		
6	PO ₄	p-Nitrophenyl-phosphate	0.1%		
7	αGLU	p-Nitrophenyl-α-D-glucoside	0.1%		
8	βGLU	p-Nitrophenyl-β-D-glucoside	0.1%	The enzymatic hydrolysis of the colorless aryl-substituted glycoside or phosphoester releases yellow σ- or p-nitrophenol, which is detected by the formation of a yellow color.	1,7,12,14, 21
9	ONPG	σ-Nitrophenyl-β-D-galactoside	0.1%		
10	GUR	p-Nitrophenyl-β-D-glucuronide	0.1%		
11	NAGA	p-Nitrophenyl-N-acetyl-β-D-glucosaminide	0.1%		
12	URE	Urea	0.6%	Utilization of urea creates basic products that change the pH indicator.	1,9,11,14,16
13	PYR	Pyrrolidine-β-naphthylamine	0.05%		
14	ARG	Arginine-β-naphthylamine	0.05%	Enzymatic hydrolysis of the arylamide substrate releases free β-naphthylamine, which is detected with RapID™ STAPH PLUS Reagent.	7,14,21,25
15	ALA	Alanine-β-naphthylamine	0.05%		
16	LEU	Leucine-β-naphthylamine	0.05%		
17	LGLY	Leucyl-glycine-β-naphthylamine	0.05%		
18	NIT	Potassium nitrate	0.015%	Utilization of nitrate results in the formation of nitrite, which is detected by nitrate reagents.	1,9,11,14,16

STORAGE

RapID™ STAPH PLUS System should be stored in its original container at 2-8°C until used. Allow product to come to room temperature before use. Remove only the number of panels necessary for testing. Reseal the plastic pouch and promptly return to 2-8°C. Panels must be used the same day they are removed from storage. RapID™ Inoculation Fluid should be stored in its original container at room temperature (20-25°C) until used.

PRODUCT DETERIORATION

This product should not be used if (1) the color of the reagent has changed, (2) the expiration date has passed, (3) the plastic tray is broken or the lid is compromised, or (4) there are other signs of deterioration.

SPECIMEN COLLECTION, STORAGE, AND TRANSPORT

Specimens should be collected and handled following recommended guidelines.^{1,12,14,17}

MATERIALS SUPPLIED

(1) 20 RapID™ STAPH PLUS System Panels, (2) 20 report forms, (3) RapID™ STAPH PLUS Reagent (one plastic dropper-bottle containing reagent sufficient for 20 panels), (4) 2 chipboard incubation trays, and (5) Instructions for Use (IFU).

MATERIALS REQUIRED BUT NOT SUPPLIED

(1) Loop sterilization device, (2) Inoculating loop, swabs, collection containers, (3) Incubators, (4) Supplemental media, (5) Quality control organisms, (6) Gram stain reagents, (7) Microscope slides, (8) Cotton swabs, (9) RapID™ Inoculation Fluid - 2 ml (REF R8325106), (10) RapID™ Nitrate A Reagent (REF R8309003), (11) RapID™ Nitrate B Reagent (REF R8309004), (12) McFarland #3 turbidity standard or equivalent (REF R20413), (13) Pipettes, and (14) ERIC® (Electronic RapID™ Compendium, REF R8323600).

PROCEDURE

NOTE: Do not interchange reagents among different RapID™ Systems.

Inoculum Preparation:

1. Test organisms must be grown in pure culture, examined by Gram stain, and tested for catalase production prior to inoculum preparation. Only catalase-positive, Gram-positive cocci that morphologically resemble staphylococci should be tested using RapID™ STAPH PLUS System. Gram-positive, catalase-negative cocci in chains, resembling streptococci, and gram-positive bacilli should not be tested.
2. Test organisms should be removed from growth media commonly used for staphylococci. The following media are recommended:
Blood Agar (TSA w/ 5% Sheep Blood), Columbia Blood Agar w/ 5% Sheep Blood, or Columbia CNA Agar, w/ Sheep Blood.

Notes:

- Plates used for inoculum preparation should be less than 72 hours old. Plates cultivated for 18-24 hours are recommended. Slow growing strains may be incubated for 48-72 hours if insufficient growth is available after 18-24 hours.
 - The use of media other than those recommended may compromise test performance.
3. Using a cotton swab or inoculating loop, suspend sufficient growth from the culture in RapID™ Inoculation Fluid (2 ml) to achieve a visual turbidity **minimally equivalent to a #3 McFarland turbidity standard.**

Suspensions less turbid than a #3 McFarland standard will result in aberrant reactions. Suspensions slightly more turbid will not compromise test performance and are recommended for stock cultures and quality control strains.

Notes:

- Suspensions should be mixed thoroughly and vortexed if required.
- Suspensions should be used within 15 minutes of preparation.


4. An agar plate may be inoculated to verify purity and any additional testing that may be required using a loopful of the test suspension from the Inoculation Fluid tube. Incubate the plate for 18-24 hours at 35-37°C.

Inoculation of RapID™ STAPH PLUS Panel:


1. Peel back the panel lid over the inoculation port by pulling the tab marked "Peel to Inoculate" up and to the left.
2. Using a pipette, gently transfer the entire contents of the Inoculation Fluid tube into the upper right-hand corner of the panel and reseal the inoculation port by pressing the peel tab back in place.

Note: It is important that the panel receive the entire test suspension.

3. After adding the test suspension and while keeping the panel on a level surface, tilt the panel back away from the test cavities at approximately a 45° angle. See below:


4. While tilted back, gently rock the panel from side to side to evenly distribute the inoculum along the rear baffles as illustrated below.

5. While maintaining a level horizontal position (best achieved by using the bench top against the reaction cavity bottoms), slowly tilt the panel forward toward the reaction cavities until the inoculum flows along the baffles into the reaction cavities (see below). This should evacuate all of the inoculum from the rear portion of the panel.

Note: If the panel is tilted too quickly, air may be trapped at the test cavity junction restricting fluid movement.

6. Return the panel to a level position. If necessary, gently tap the panel on the bench top to remove any air trapped in the cavities.

Notes:

- Examine the test cavities to verify they are bubble-free and uniformly filled. Slight irregularities in test cavity fill are acceptable and will not affect test performance. If the panel is grossly misfilled, a new panel should be inoculated and the misfilled panel discarded.
- Complete the inoculation of each panel receiving inoculation fluid before inoculating additional panels.
- Do not let the inoculum rest in the back portion of the panel for prolonged periods without completing the procedure.
- Used Inoculation Fluid tubes, swabs, and other contaminated materials should be properly sterilized prior to disposal.

Incubation of RapID™ STAPH PLUS Panel:

Incubate inoculated panels at 35-37°C in a non-CO₂ incubator for **at least 4 HOURS and not more than 6 HOURS**. For ease of handling, panels may be incubated in the chipboard incubation trays provided with the kit.

Scoring of RapID™ STAPH PLUS Panels:

RapID™ STAPH PLUS Panels contain 18 reaction cavities that provide 18 test scores. Reagent-requiring tests (cavities 13-18) are designated with a box drawn around them as illustrated below. Refer to Table 2 for complete information on interpreting RapID™ STAPH PLUS System tests.

RapID™ STAPH PLUS Panel Test Location:

Cavity	1	2	3	4	5	6	7	8	9	10	11	12	13			14		15		16		17		18	
Test Code	ADH	ODC	LIP	SUC	MANO	PO ₄	αGLU	βGLU	ONPG	GUR	NAGA	URE	PYR			ARG		ALA		LEU		LGLY		NIT	
Reagent	None												Rapid™ STAPH PLUS Reagent						Rapid™ Nitrate A Reagent		Rapid™ Nitrate B Reagent				

- While firmly holding the RapID™ STAPH PLUS Panel on the bench top, peel off the label lid over the reaction cavities by pulling the lower right hand tab up and to the left.
 - Cavities 1 (ADH) through 12 (URE) are read without addition of reagent. To the remaining cavities (13 through 18), add the following reagents as indicated:
 - Add 2 drops of RapID™ STAPH PLUS Reagent to cavities 13 (PYR) through 17 (LGLY).
 - Add 1 drop of RapID™ Nitrate A Reagent to cavity 18 (NIT).
 - Add 1 drop of RapID™ Nitrate B Reagent to cavity 18 (NIT).
- NOTE:** In most cases, immediate color development will be observed after addition of reagents. If immediate color development is not observed, allow approximately 30 seconds (but no more than 1 minute) for color development.
- After reagent addition, read and score the test cavities from left to right using the interpretation guide presented in Table 2. Record the scores in the appropriate boxes on the report form.
 - To identify the isolate, reference the microcode obtained in ERIC® or use the Differential Chart included in this IFU.

Table 2. Interpretation of RapID™ STAPH PLUS Panel Tests*

Cavity #	Reagent	Test Code	Reaction		Comment
			Positive	Negative	
1 2	None	ADH ODC	Purple, blue, bluish gray, or greenish gray	Yellow or, straw	Any purple or significant blue, or gray color should be scored positive.
3 4 5	None	LIP SUC MANO	Yellow, yellow orange	Red, dark red orange, or orange	Any yellow or yellow orange throughout the well should be scored positive.
6	None	PO ₄	Yellow or pale yellow	Clear, cream, or very light yellow	Only significant yellow color throughout the well should be scored positive. Very light yellow should be scored NEGATIVE.
7 8 9 10 11	None	αGLU βGLU ONPG GUR NAGA	Any yellow or light yellow	Clear, cream or buff	Any yellow color throughout the well should be scored positive. Note, yellow color may be faint with some taxa.
12	None	URE	Red or dark red orange	Yellow, yellow orange, or very pale orange	Any red or dark orange throughout the well should be scored positive.
13	RapID™ STAPH Plus Reagent	**PYR	Purple or dark red	Yellow, orange, pale light red or shades of pink	Only significant purple or red should be scored positive.
14		ARG	Purple, red, light red or dark pink	Yellow, orange or very pale pink	Any purple, red or dark pink should be scored positive. Yellow, pale light pinks, and orange are scored negative.
15		ALA			
16		LEU			
17		LGLY			
18	RapID™ Nitrate A and Nitrate B	NIT	Cherry red or dark pink	Clear, yellow, cream or very pale pink	Only a significant red or dark pink, developing immediately or within 30 seconds, should be scored positive. Pinks that develop after prolonged periods after reagent addition should be scored negative.

*NOTE: Panels are best observed by looking down through the reaction cavities against a white background.

** NOTE: PYR color development is interpreted differently from other tests requiring addition of RapID™ STAPH PLUS Reagent.

RESULTS AND RANGE OF EXPECTED VALUES

The RapID™ STAPH PLUS Differential Chart illustrates expected results for the RapID™ STAPH PLUS System. Differential Chart results are expressed as a series of positive percentages for each system test. This information statistically supports the use of each test and provides the basis, through numerical coding of digital test results, for a probabilistic approach to the identification of the test isolate.

Identifications are made using individual test scores from RapID™ STAPH PLUS panels in conjunction with other laboratory information (e.g., Gram stain, colonial morphology, coagulase test) to produce a pattern that statistically resembles known reactivity for taxa recorded in the RapID™ STAPH PLUS System database. These patterns are compared through the use of the RapID™ STAPH PLUS Differential Chart or by derivation of a microcode and the use of ERIC®.

LIMITATIONS

1. The use of RapID™ STAPH PLUS System and the interpretation of results require a competent laboratorian who is trained in general microbiological methods and who judiciously makes use of knowledge, experience, specimen information, and other pertinent procedures before reporting the isolate identity obtained.
2. The RapID™ STAPH PLUS System must be used with pure cultures of test organisms. The use of mixed microbial populations or direct testing of clinical material without culture will result in aberrant results.
3. The RapID™ STAPH PLUS System is designed for use with the taxa listed in the RapID™ STAPH PLUS Differential Chart. Gram-positive bacilli, gram-negative cocci, or catalase-negative gram-positive cocci in chains should not be tested.
4. Expected values listed for RapID™ STAPH PLUS System tests may differ from conventional test results or previously reported information.
5. The accuracy of the RapID™ STAPH PLUS System is based upon the statistical use of a multiplicity of specially designed tests and an exclusive, proprietary database. The use of any single test found in the RapID™ STAPH PLUS System to establish the identification of a test isolate is subject to the error inherent in that test alone.
6. Results obtained using the RapID™ STAPH PLUS System are dependent upon adherence to the procedures indicated. Any change or modification in the procedure may result in aberrant results.

PERFORMANCE CHARACTERISTICS

The RapID™ STAPH PLUS System performance characteristics have been established by laboratory testing of clinical, reference, and stock cultures at Remel. Among 291 strains tested, 280 (96%) agreed with the reported reference identifications (with or without supplementary tests). One (0.3%) strain provided a questionable microcode that did not result in an identification. Ten strains (3.4%) did not agree with reported reference identifications.

BIBLIOGRAPHY

1. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. *Bailey and Scott's Diagnostic Microbiology*. 9th ed. Mosby Year Book, Inc., St. Louis, MO.
2. Buttery, J.P., M. Easton, S.R. Pearson, and G.G. Hogg. 1997. *J. Clin. Microbiol.* 35:2174-2177.
3. Ben-Ami, R., S. Navon-Venezia, D. Schwartz, and Y. Carmeli. 2003. *J. Clin. Microbiol.* 41:2444-2447.
4. Calvo, J., J.L. Hernández, M.C. Fariñas, D. García-Palomo, and J. Agüero. 2000. *J. Clin. Microbiol.* 38: 3887-3889.
5. De Paulis, A.N., S.C. Predari, C.D. Chazarreta, and J.E. Santoianni. 2003. *J. Clin. Microbiol.* 41:1219-1224.
6. Funke, G. and P. Funke-Kissling. 2005. *J. Clin. Microbiol.* 43:84-88.
7. Guilbert, G.G. 1970. *Methods of Enzymatic Analysis*. Pergamon Press, New York, NY.
8. Haile, D.T., J. Hughes, E. Vetter, P. Kohner, R. Snyder, R. Patel, and F. R. Cockerill III. 2002. *J. Clin. Microbiol.* 40:654-656.
9. Heikens, E., A. Fleer, A. Paauw, A. Florijn, and A.C. Fluit. 2005. *J. Clin. Microbiol.* 43:2286-2290.
10. Horstkotte, M.A., J.K. Knobloch, H. Rohde, S. Dobinsky, and D. Mack. 2004. *J. Clin. Microbiol.* 42:5041-5046.

11. Ieven, M., J. Verhoeven, S.R. Pattyn, and H. Goossens. 1995. *J. Clin. Microbiol.* 33:1060-1063.
12. Isenberg, H.D., 1992. *Clinical Microbiology Procedures Handbook*. Vol.1. ASM, Washington, D.C.
13. Kawamura, Y., X.G. Hou, F. Sultana, K. Hirose, M. Miyake, S. E. Shu, and T. Ezaki. 1998. *J. Clin. Microbiol.* 36:2038-2042.
14. Murray, P.R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R.H. Tenover. 1999. *Manual of Clinical Microbiology*. 7th ed. ASM, Washington, D.C.
15. Mahoudeau, I., X. Delabranche, G. Prevost, H. Monteil, and Y. Piemont. 1997. *J. Clin. Microbiol.* 35:2153-2154.
16. MacFaddin, J.F. 1980. *Biochemical Tests for Identification of Medical Bacteria*. 2nd ed. Williams and Wilkins, Baltimore, MD.
17. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenover. 2003. *Manual of Clinical Microbiology*. 8th ed. ASM, Washington, D.C.
18. Pottumarthy, S., J.M. Schapiro, J.L. Prentice, Y.B. Houze, S.R. Swanzey, F.C. Fang, and B.T. Cookson. 2004. *J. Clin. Microbiol.* 42: 5881-5884.
19. Renneberg, J., K. Rieneck, and E. Gutschik. 1995. *J. Clin. Microbiol.* 33:1150-1153.
20. Rhoden, D.L. and J.M. Miller. 1995. *J. Clin. Microbiol.* 33:96-98.
21. Bascomb, S. and M. Manafi. 1998. *Clin. Microbiol. Rev.* 11:318-340.
22. Shuttleworth, R., R.J. Behme, A. McNabb, and W.D. Colby. 1997. *J. Clin. Microbiol.* 35:2537-2541.
23. Srinivasan, A., J.D. Dick, and T.M. Perl. 2002. *Clin. Microbiol. Rev.* 15:430-438.
24. Stepanovic, S., I. Dakic, D. Morrison, T. Hauschild, P. Jezek, P. Petras, A. Martel, D. Vukovic, A. Shittu, and L.A. Devriese. 2005. *J. Clin. Microbiol.* 43:956-958.
25. Norris, J.R. and D.W. Ribbons. 1976. *Methods in Microbiology*. Vol. 9. Academic Press, New York, NY.
26. Weinstein, M.P., S. Mirrett, L. Van Pelt, M. McKinnon, B. L. Zimmer, W. Kloos, and L.B. Reller. 1998. *J. Clin. Microbiol.* 36:2089-2092.

PACKAGING


REF R8311009, RapID™ STAPH PLUS System20 Tests/Kit

Symbol Legend

REF	Catalog Number
IVD	In Vitro Diagnostic Medical Device
LAB	For Laboratory Use
	Consult Instructions for Use (IFU)
	Temperature Limitation (Storage Temp.)
LOT	Batch Code (Lot Number)
	Use By (Expiration Date)
EC REP	European Authorized Representative

RapID™ is a trademark of Remel Inc.
 ERIC® is a registered trademark of Remel Inc.
 ATCC® is a registered trademark of American Type Culture Collection.

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RapID™ es una marca comercial de Remel Inc.
ERIC® es una marca registrada de Remel Inc.