# mEC Broth with Novobiocin

For the selective enrichment of enterohemorrhagic E. coli (EHEC) in foods.

The culture medium complies with the requirements of the USDA-FSIS method for the isolation and identification of Escherichia coli 0157:H7 from meat.

# Mode of Action

The nutrient substrates contained in mEC Broth provide favourable growth conditions. Especially lactose improves proliferation of lactose-positive bacteria. The mixture of bile salts No. 3 and Novobiocin suppresses the growth of Grampositive microbial flora.

# Typical Composition (g/litre)

Peptone 20.0; sodium chloride 5.0; bile salts No. 3 1.12; lactose 5.0; di-potassium hydrogen phosphate 4.0; potassium dihydrogen phosphate 1.5; novobiocin 0.02.

# Preparation

Suspend 36.7 g in 1 litre of demin. water; autoclave (15 min at 121  $^{\circ}$ C).

pH: 6.9 ± 0.2 at 25 °C.

The prepared broth is clear and yellow-orange. The broth is stable for up to 6 months when stored at +2 to +8  $^\circ$ C.

#### Experimental Procedure and Evaluation

Inoculate mEC Broth with the sample material (usually add 25 g test portion of the sample to 225 ml of broth).

Incubation: 18 to 24 h at 37 °C or 41.5 °C. Incubation temperature depends on standard used. Afterwards about 0.1 ml of the broth is streaked on the dry surface of a E. coli 0157:H7 Selective Agar, e.g. Fluorocult® E. coli 0157:H7 Agar, SMAC Agar or CT-SMAC Agar in a way that single colonies can be well isolated.

### Quality control

Test strains	Growth	Singlepath® E. coli 0157
E. coli 0157:H7 ATCC 35150 (0157)	good	+
E. coli ATCC 11775	fair / good	-
Proteus vulgaris ATCC 13315	good inhibition	
Pseudomonas aeruginosa ATCC 27853	good inhibition	
Enterococcus faecalis ATCC 33186	good inhibition	

# Literature

OKREND, A.J.G., ROSE, B.E., a. BENNETT, B.: A research not: A screening method for the isolation of Escherichia coli 0157:H7 from ground beef. -J.Food Prot., 53; 249-252 (1990).

USDA-FSIS, Revision 4 of Laboratory Communication #38 Protocol for Isolation and Identification of Escherichia coli 0157:H7. - Amelia K. Sharar and Bonnie E. Rose, (1996).

#### **Ordering Information**

Product	Merck Cat. No.	Pack size
mEC Broth with Novobiocin	1.14582.0500	500 g
CT-Supplement	1.09202.0001	1 x 16 vials
Fluorocult® E. coli 0157:H7 Agar	1.04026.0500	500 g
Sorbitol-MacConkey Agar	1.09207.0500	500 g
Singlepath <sup>®</sup> E.coli 0157	1.04141.0001	25 tests