

# Combi K-Set


**Rapid diagnostic test for *in vitro* detection of Rotavirus & Adenovirus in stool specimens**

- Safe
- Easy to use
- Ready to use with a complete set of sampling devices if needed
- Accurate: same performances as dipsticks
- Easy to interpret



[www.corisbio.com](http://www.corisbio.com)

Accurate, economical and fast procedure!



## Performance compared to EIA

	Rotavirus detection (n=214)	Adenovirus detection (n=130)
Specificity	100% (108/108)	100% (118/118)
Sensitivity	99.1% (105/106)	100% (12/12)
PPV	100% (105/105)	100% (12/12)
NPV	99.1% (108/109)	100% (118/118)
Accuracy	99.5% (213/214)	100% (130/130)

## Bibliography

For more information:

[www.corisbio.com](http://www.corisbio.com) or [info@corisbio.com](mailto:info@corisbio.com)

## Ordering Information

Immunochromatographic range

### Enteric (20 tests)

Description	Code
Combi K-SeT	K-1504*
Giardia K-SeT	K-1513*

### Gastric (20 tests)


Description	Code
Pylori K-SeT	K-1519*

### Respiratory (20 tests)

Description	Code
RSV Respi K-SeT	K-1506*
Adeno Respi K-SeT	K-1509*
Influ-A Respi K-SeT	K-1510*

\* Also available in dipstick

Distributed by:



Science Park CREALYS  
Rue Jean Sonet 4A - 5032 Gembloux - Belgium  
Ph.: +32 (0)81 719 917 - Fx.: +32 (0)81 719 919  
[info@corisbio.com](mailto:info@corisbio.com) - [www.corisbio.com](http://www.corisbio.com)

MADE IN BELGIUM

Fo 7.5-34-04/00  
October 2010